Bilirubin Control – Pediatric / Level 1 & 2

IVD 131651 131652 REF 1331-31 1332-31 2023-03-31 CE IVD 2°C 8°C i

CE IVD IVD Cont. REF i EC REP i i

European Conformity CE-Konformitätskennzeichnung Conformité aux normes européennes Conformité europea Conformidad europea
Lot Number Bezeichnung Designation du lot Conformité europea Denominación de lote
Manufactured by Hersteller von Fabriqué par Fabricado por
For in vitro diagnostic use In-vitro Diagnosticum Pour diagnostic in vitro Per uso diagnostico in vitro De uso diagnostico in vitro
Biological Risk Biogefährdung Risque biologique Riesgo biológico
Contents of kit Inhalt der Packung Contenu du coffret Contenido della confezione Contenido del estuche
Catalog No. Bestellnr. N° de catalogue Catalogo n. Nº de catálogo
Caution, See Product Insert Achtung, Siehe Packungsbeilage Attention, voir notice d'utilisation Attenzione, vedere il foglio illustrativo del prodotto Atención, consulte el folleto del producto
Authorized Representative Bevollmächtigter Représentant agréé Representante autorizzato Representante autorizado
Temperature Limitation Temperaturbegrenzungen Limites de température Limiti di temperatura Limite de temperatura
Consult instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation Consultare le istruzioni l'uso Consultare le istruzioni de uso
Use by (last day of month) Verwendbar bis (letzter Tag des Monats) Utilisable jusqu'à (dernier jour du mois inclus) Da utilizzare prima del (ultimo giorno del mese) Estable hasta (ultimo día del mes)

English Intended Use

The Quantimetrix Pediatric Bilirubin Control is intended as a means of monitoring various bilirubin assay methods to validate quantitation of patient samples.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Quantimetrix Pediatric Bilirubin Control is supplied as a liquid, 3 x 3 mL vials per box. It is prepared from purified bilirubin in a human protein base. Stabilizers and preservatives have been added.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the protein pool have been tested and found non-reactive for Hepatitis B surface antigen, and for HCV antibody and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is suggested that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Storage and Stability

The control should be stored at 2°C–8°C. When stored unopened at 2°C–8°C, the control is stable until the expiration date stated on the label. If the rubber stopper is removed or air is introduced into the vial, the stability of the control is limited to twenty-one (21) days at 2°C–8°C. **Longer stability may be achieved by removing the sample with a syringe.** Discard the control if it is turbid, if there is any evidence of microbial contamination or if it exhibits any green color. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the control from the refrigerator and allow it to come to room temperature (18°C–25°C) for 10 to 15 minutes. Gently invert the control to assure homogeneity of the contents. Avoid foaming the control. For maximum stability, lift up the center tab of the metal seal without removing the seal. Pierce the rubber stopper with a syringe and withdraw the required quantity of control. **Do not introduce air into the vial.** Do not return unused material to the vial. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method.

Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

Limitations

Bilirubin is sensitive to light and oxygen. Exposure to these factors may decrease the values from those stated in the expected ranges.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

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Analytes / Instruments	Units	Level 1 / Lot 131651		Level 2 / Lot 131652		SI Units ¹	Level 1 / Lot 131651		Level 2 / Lot 131652	
		mean	expected range	mean	expected range		mean	expected range	mean	expected range
Total Bilirubin										
Abbott Architect*	mg/dL	10.5	8.4 – 12.6	20.1	16.1 – 24.1	µmol/L	180	144 – 216	344	275 – 413
Siemens Dimension	mg/dL	10.5	8.4 – 12.6	20.2	16.2 – 24.3	µmol/L	179	143 – 215	346	277 – 415
Reichert UNISTAT® Bilirubinometer	mg/dL	11.6	9.3 – 13.9	23.8	19.1 – 28.6	µmol/L	198	158 – 237	408	326 – 489
Beckman Coulter AU400/AU640/AU2700*	mg/dL	9.1	7.2 – 10.9	18.4	14.7 – 22.1	µmol/L	155	124 – 186	315	252 – 377
Advanced Instruments	mg/dL	10.8	8.7 – 13.0	21.3	17.1 – 25.6	µmol/L	185	148 – 222	365	292 – 438
Roche Hitachi (Roche Reagent)	mg/dL	10.5	8.4 – 12.6	20.7	16.5 – 24.8	µmol/L	180	144 – 216	353	282 – 424
Direct Bilirubin										
Abbott Architect*	mg/dL	2.4	1.9 – 2.9	4.7	3.7 – 5.6	µmol/L	41	33 – 49	80	64 – 96
Siemens Dimension	mg/dL	2.3	1.8 – 2.8	4.6	3.7 – 5.5	µmol/L	39	31 – 47	78	63 – 94
Beckman Coulter AU400/AU640/AU2700*	mg/dL	6.5	5.2 – 7.7	9.4	7.5 – 11.2	µmol/L	110	88 – 132	160	128 – 192
Advanced Instruments	mg/dL	2.8	2.2 – 3.4	5.3	4.3 – 6.4	µmol/L	48	38 – 57	91	73 – 110
Roche Hitachi (Roche Reagent)	mg/dL	2.8	2.2 – 3.3	5.6	4.4 – 6.7	µmol/L	48	38 – 57	95	76 – 114
Neonatal										
DAS Neo Bil Plus	mg/dL	10.9	8.7 – 13.1	20.3	16.3 – 24.4	µmol/L	186	149 – 223	348	278 – 417

Footnotes

¹ SI International System of Units
* Mean and range based on limited data

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Abbott Laboratories, Abbott ARCHITECT®
Advanced Instruments, BR2 Bilirubin STAT-Analyzer™
Beckman Coulter, Inc., AU Clinical Chemistry System

Reichert UNISTAT® Bilirubinometer
Roche Diagnostics Corporation, Roche Cobas® 6000, COBAS® INTEGRA 400 Plus
Siemens Healthcare, Dimension® System, Dimension® EXL™ & Dimension Vista® Systems

Quantimetrix Corporation
2005 Manhattan Beach Blvd.
Redondo Beach, CA 90278
+1.310.536.0006

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MDSS GmbH