
# Quantimetrix® CHROMASCOPICS®

## Urinalysis Control with Microscopics Level 1 & 2

**LOT** 237141  
237142 **REF** 1540-01  
1540-03 **EXP** 2025-11-30

### English

#### Intended Use

The Quantimetrix **CHROMASCOPICS** Urinalysis Control with Microscopics is intended as an assayed quality control for the Siemens Clinitek Novus Automated Urine Chemistry Analyzer and Atellica 800 and 1500 systems. It is also compatible with Siemens urinalysis dipsticks, confirmatory tests such as **K-CHECK** and **Ictotest**® reagent tablets, and qualitative **hCG** methods.

In addition, **CHROMASCOPICS** may also be used as a means of validating the processing and centrifugation of patient urine samples prior to the microscopic evaluation of urine sediment. For professional use only.

#### Summary and Explanation

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

Microscopic QC controls must be run each day the test is performed.<sup>1</sup> Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.<sup>2</sup> Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts, and crystals.<sup>3,4</sup>

#### Product Description

**CHROMASCOPICS** is supplied as a liquid, ready-to-use control in two levels, configured with 2x120 mL of each level per box or 2x25 mL per bag. It does not require reconstitution or dilution. They are prepared from human urine to which stabilized human red and white blood cells, calcium oxalate crystals, non-pathogenic bacteria, and other compounds have been added to produce the desired reactions when tested by the methods indicated in the **Intended Use** section. Preservatives have been added to inhibit microbial growth.

#### Caution

POTENTIAL BIOHAZARDOUS MATERIAL. Contains human urine, human blood cells and human Chorionic Gonadotropin (hCG) from pregnancy urine. The human hCG source material and all blood donor units comprising the human cell source material used in the manufacture of this product have been tested and found nonreactive for Hepatitis B Surface Antigen and Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

#### Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, 1,3,3-Trimethyl-2-[(2-methyl-3H-indol-3-ylidene)ethylidene]indoline monohydrochloride, 1,2-Propylene Glycol, Level 1; 2,4-Pentanedione, Level 2.

H317 – May cause an allergic skin reaction.

P261 – Avoid breathing vapors, mist, or spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations. Safety Data Sheet (SDS) available for professional users at [quantimetrix.com](http://quantimetrix.com).

#### Storage and Stability

The **CHROMASCOPICS** Urinalysis Control with Microscopics Kit should be stored at 2°C–8°C when not in use. **Do not freeze.** When stored at 2°C–8°C the controls are stable until the expiration date stated on the label. After opening, the controls will remain stable until the expiration date stated on the label when stored at 2°C–8°C between uses. Discard the control if it becomes more turbid or develops a stronger odor. Discard controls in the same manner as other biological specimens, according to local guidelines.

Remove the controls from the refrigerator and replace the cap on the control bottle with the spout cap included in the kit. Allow the control to come to room temperature (18°C–25°C) for approximately 15–30 minutes depending on the volume remaining in the bottle. Mix the controls thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Immediately after dispensing close the spout cap and store the controls at 2°C–8°C when not in use.

#### Procedure for Urinalysis and Microscopic Evaluation of Urine Sediment

For the Siemens Clinitek Novus, Atlas, and Atellica 800 and 1500 systems follow the manufacturer's instructions for control material. For Siemens Semi-Automated and POC urinalysis systems and for visual testing dispense the control into a tube and test as if it was a patient specimen. Read the urinalysis reagent strips, visually or with an instrumental reader, in accordance with the manufacturer's instructions.

**For manual microscopic evaluation of urine sediment**, treat the controls as you would patient samples in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.<sup>1</sup>

#### Procedure for hCG Tests and Confirmatory Tests

**Note:** The bottles of Level 1 Control are to be used as negative controls for hCG methods. The bottles of Level 2 Control are to be used as positive controls for hCG methods. Use the negative and positive controls as if they were patient specimens in accordance with the test kit manufacturer's instructions.

Most manufacturers of pregnancy test kits specify the volume of sample to be used with their kits. Many kits include transfer pipets to be used to deliver a certain sample volume onto the test device. It is important that sufficient volume be used to produce the correct test result. If dispensing the control for hCG tests and confirmatory tests directly from the control bottle, each user should validate that the volume (number of drops) dispensed by the included spout cap is sufficient to meet the pregnancy test kit's and confirmatory tests' requirement for sample volume.

If using the same bottle of control dispensed for urinalysis testing and microscopic evaluation, remove the volume of sample to be used for hCG tests and confirmatory tests after centrifugation, before discarding the supernatant and without disturbing the sediment.

#### Expected Values

All expected value ranges and assignments are determined using pooled data from the system manufacturer and/or internal data generated as per the manufacturer's instructions.

#### Limitations

After running **CHROMASCOPICS** on the Siemens Clinitek Novus Automated Urine Chemistry Analyzer, the operator must clean the SG well. This is assessable via the System Function Bar on the home screen. Refer to the Siemens Clinitek Novus Automated Urine Chemistry Analyzer Operators Manual for detailed instructions.

**CHROMASCOPICS** contains a non-hazardous dye. If spilled it may stain some surfaces. Bleach will remove the dye from most surfaces.

**CHROMASCOPICS** is primarily intended for use with the Siemens Clinitek Novus Automated Urine Chemistry Analyzer and Atellica 800 and 1500 systems. It also has been demonstrated to be compatible with the other products as stated in the **Intended Use** section. For all other urinalysis test methods please consult the product descriptions on our web site.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Detailed information on the limitations of each test method is included in the limitations section of the manufacturers' package insert. Technical updates can be found on our website. The Quality Control Log can be downloaded from the Quantimetrix website at [quantimetrix.com](http://quantimetrix.com) or contact Tech Support at (310) 536-0006, option 3.

The appearance of a macroscopic crystalline precipitate in the product will not affect performance.

#### Multistix Users

Colors produced by the **urobilinogen** and/or **bilirubin** reactions on these dipsticks with the Urinalysis Dipstick Control may not be characteristic of those shown on the manufacturer's label when reading the dipstick reactions visually. The urobilinogen reactions are consistent and intensify with the increase in the urobilinogen concentration but may not provide an exact color match to those displayed on the label.

**Note:** Siemens® CLINITEK 50 and Siemens® STATUS or CLINITEK STATUS PLUS may see an Albumin/ Creatinine ratio result of "Abnormal" with the Level 1 control.

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- 1 Routine Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, NCCLS Documentation GP 16-A, Approved Guideline, (1995), NCCLS, Wayne PA, 19087.
- 2 Bradley GM, Benson ES, Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed., Philadelphia, PA: Saunders, 1974.
- 3 Bologna, CV, Understanding Laboratory Medicine, St. Louis, MO: CV Mosby Co., 1971.
- 4 Dudas HC, Lab Med, 12:765, 1981.

Analytes/Method	Level 1 - 237141	Level 2 - 237142	Units
<b>Red Blood Cells (Erythrocytes)</b>			
Siemens Atellica 800/1500 Systems	3 - 17	36 - 90	Cells/ $\mu$ L
Siemens Atellica 800/1500 Systems	1 - 4	8 - 20	Cells/hpf
KOVA® Glasstic® Slide 10 with Grids	0 - 16	40 - 123	p/ $\mu$ L
Non-grid slides (~0.5 mL)	0 - 6	5 - 54	p/hpf
Non-grid slides (~1.0 mL)	0 - 4	4 - 25	p/hpf
FisherBrand UriSystem DeciSlide	0 - 5	3 - 21	p/hpf
<b>White Blood Cells (Leukocytes)</b>			
Siemens Atellica 800/1500 Systems	2 - 23	23 - 73	Cells/ $\mu$ L
Siemens Atellica 800/1500 Systems	0 - 5	5 - 16	Cells/hpf
KOVA® Glasstic® Slide 10 with Grids	0 - 10	9 - 62	p/ $\mu$ L
Non-grid slides (~0.5 mL)	0 - 7	2 - 22	p/hpf
Non-grid slides (~1.0 mL)	0 - 4	1 - 14	p/hpf
FisherBrand UriSystem DeciSlide	0 - 4	1 - 12	p/hpf

Analytes/Method	Level 1 - 237141	Level 2 - 237142
<b>Casts</b>		
Siemens Atellica 800/1500 Systems	none	none
KOVA® Glasstic® Slide 10 with Grids	none	none
Non-grid slides (~0.5 mL)	none	none
Non-grid slides (~1.0 mL)	none	none
FisherBrand UriSystem DeciSlide	none	none
<b>Crystals (calcium oxalate dihydrate with some monohydrate and amorphous forms)<sup>12</sup></b>		
Siemens Atellica 800/1500 Systems	negative	1+ - 3+
KOVA® Glasstic® Slide 10 with Grids	none	present
Non-grid slides (~0.5 mL)	none	present
Non-grid slides (~1.0 mL)	none	present
FisherBrand UriSystem DeciSlide	none	present
<b>Bacteria</b>		
Siemens Atellica 800/1500 Systems	0 - 10	934 - 2268
KOVA® Glasstic® Slide 10 with Grids	may be present	present
Non-grid slides (~0.5 mL)	may be present	present
Non-grid slides (~1.0 mL)	may be present	present
FisherBrand UriSystem DeciSlide	may be present	present

Analytes	Level 1 - 237141	Level 2 - 237142
<b>Confirmatory and Other Tests</b>		
K-CHECK (Ketones)	Negative	Small - Large
Ictotest (Bilirubin)	Negative	Positive
Refractometer (Specific Gravity)	1.019 - 1.025	1.011 - 1.017
hCG	Negative	Positive
pH Paper	4.0 - 6.0	7 - 9
Sulfosalicylic Acid (Total Protein)	Negative ( $\leq 0.05$ ) <sup>10</sup>	Positive ( $\geq 0.50$ ) <sup>10</sup>

Analytes	Level 1 - 237141	Level 2 - 237142
<b>SIEMENS VISUAL TESTING (Visual Test Strips Only)</b>		
Glucose	Negative	100 - 1000 mg/dL
Bilirubin	Negative	Sm - Lg (1+ - 3+)
Ketones	Negative	Trace (5) - 160 mg/dL (Tr - Lg)
Specific Gravity	1.015 - 1.025	1.000 - 1.015
Blood	Negative - Trace	Sm - Lg (1+ - 3+)
pH	5 - 6.5	7.5 - $\geq 9.0$
Protein	Negative	Trace - $\geq 2000$ mg/dL (Tr - 4+)
Urobilinogen	Normal (0.2 E.U./dL) <sup>6</sup>	2.0 - 8.0 E.U./dL *
Nitrite	Negative	Positive
Leukocytes	Negative	Tr - Lg (Tr - 3+)
Creatinine <sup>3</sup>	10 - 50 mg/dL	100 - 300 mg/dL
<b>SIEMENS CLINITEK 50</b>		
Glucose	Negative	100 - 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - $\geq 80$ mg/dL ( $\pm$ - 3+)
Specific Gravity	1.010 - $\geq 1.030$	1.000 - 1.020
Blood	Negative - Trace	Trace - Large (Tr - 3+)
pH	5.0 - 6.0	7.0 - 8.5
Protein	Negative	30 - $\geq 300$ mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - 4.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Microalbumin <sup>1</sup>	10 - 30 mg/L	80 - 150 mg/L
Creatinine <sup>3</sup>	10 - 100 mg/dL	100 - 300 mg/dL
<b>SIEMENS CLINITEK 500</b>		
Glucose	Negative	100 - 500 mg/dL (Tr - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - $\geq 80$ mg/dL ( $\pm$ - 3+)
Specific Gravity	1.015 - $\geq 1.030$	$\leq 1.005$ - 1.020
Blood	Negative - Trace	Small - Large (1+ - 3+)
pH	5.0 - 6.5	7.5 - $\geq 9.0$
Protein	Negative	30 - $\geq 300$ mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	2.0 - 6.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Creatinine <sup>2</sup>	10 - 100 mg/dL	50 - 300 mg/dL
<b>SIEMENS CLINITEK ADVANTUS</b>		
Glucose	Negative	100 - $\geq 1000$ mg/dL (Tr - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Tr - $\geq 80$ mg/dL (Tr - 3+)
Specific Gravity	1.015 - $\geq 1.030$	$\leq 1.005$ - 1.025
Blood	Negative - Trace	Small - Large (1+ - 3+)
pH	5.0 - 6.5	7.5 - $\geq 9.0$
Protein	Negative	30 - $\geq 300$ mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	2.0 - $\geq 8.0$ E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Creatinine <sup>2</sup>	10 - 50 mg/dL	50 - 300 mg/dL
<b>SIEMENS CLINITEK STATUS / STATUS PLUS / STATUS CONNECT</b>		
Glucose	Negative	100 - 500 mg/dL (Tr - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Tr - $\geq 160$ mg/dL (Tr - 4+)
Specific Gravity	1.010 - $\geq 1.030$	1.010 - 1.025
Blood	Negative - Trace	Small - Large (1+ - 3+)
pH	5.0 - 6.0	7.0 - $\geq 9.0$
Protein	Negative	Trace - $\geq 300$ mg/dL (Tr - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	2.0 - $\geq 8.0$ E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Microalbumin <sup>1</sup>	10 - 30 mg/L	80 - 150 mg/L
Creatinine <sup>3</sup>	10 - 100 mg/dL	100 - 300 mg/dL
hCG	Negative	Positive
<b>SIEMENS CLINITEK NOVUS/ATELLICA 1500 SYSTEM</b>		
Albumin	10 - 80 mg/L	80 - 150 mg/dL
Bilirubin	Negative	Small - Large (1+ - 3+)
Blood	Negative	Small - Large (1+ - 3+)
Creatinine	10 - 50 mg/dL	50 - $\geq 300$ mg/dL
Glucose	Negative	100 - $\geq 1000$ mg/dL (Tr - 3+)
Ketones	Negative	Tr - $\geq 160$ mg/dL (Tr - 4+)
Leukocytes	Negative	Trace - Large (Tr - 3+)
Nitrite	Negative	Positive
pH	5.0 - 6.0	7.5 - $\geq 9.0$
Protein	Negative	30 - 300 mg/dL (1+ - 3+) <sup>6</sup>
Specific Gravity	1.019 - 1.025	1.012 - 1.018
Urobilinogen	Normal (0.2 E.U./dL)	0.2 - 4.0 E.U./dL
Color	$\leq$ Dark Yellow	$\leq$ Orange
Clarity	Clear - Slightly Cloudy	Clear - Slightly Cloudy

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**ENGLISH**

- 2 Values only apply to Clinitek Microalbumin Reagent Strips when read on the Clinitek 50 and Status.
- 3 Values only apply to Multistix Pro™ Reagent Strips
- 4 Values only apply to Multistix Pro and Clinitek Microalbumin Reagent Strips when read on Clinitek Urine Analyzers
- 6 Some customers may obtain false negatives.
- 7 Some customers may obtain false positives.
- 10 Absorbance at 620 nm
- 12 Calcium oxalate dihydrate with some monohydrate and amorphous forms
- \* See Limitations

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