Complete D® 25-OH Vitamin D Control / Level 1 & 2

LOT 169161
169162

REF 1290-01

2025-02-28

2°C - 8°C

European Conformity
CE-Konformitätszeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Beschreibung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabriqué par
Fabricado da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
Biogefährdung
Risque biologique
Rischio biologico
Perigo biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Besteller
N° de catalogue
Catalogo n.
N° de catálogo

Caution, See Product Insert
Achtung, Siehe Packungsbeilage
Attention, voir notice d'utilisation
Attenzione, vedere il foglietto
illustrativo del prodotto
Atención, consulte el folleto
del producto

Authorized Representative
Bevollmächtigter
Représentant agréé
Rappresentante autorizzato
Representante autorizado

Temperature limitation
Temperaturbegrenzung
Limites de température
limiti di temperatura
límite de temperatura

Consult instructions for use
Gebrauchsanweisung beachten
Consulter les instructions d'utilisation
Consultare le istruzioni l'uso
Consulte las instrucciones de uso

Use by (last day of month)
Verwendbar bis (letzter Tag des Monats)
Utilisable jusqu'à (dernier jour du mois inclus)
Da utilizzare prima del (ultimo giorno del mese)
Estable hasta (último día del mes)

English

Intended Use

The Complete D 25-OH Vitamin D Control is intended to monitor the performance of clinical assays used to quantitate total 25-OH Vitamin D.

Use of quality control materials is an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Complete D 25-OH Vitamin D Control is supplied in two levels, 3 mL per vial. Each kit contains 2 vials of level 1 and 2 vials of level 2 controls. The controls are liquid, ready to use and require no reconstitution or dilution. The controls are prepared in a human serum matrix fortified with analytes and preservatives. This control contains both 25-OH Vitamin D2 and D3; and when combined represent clinically relevant levels of Total Vitamin D.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. The Complete D Controls are prepared in a human serum matrix and contain other human source materials. All blood donor units comprising the serum pool have been tested and found nonreactive for Hepatitis B surface antigen, Hepatitis C and HIV 1 & 2 antibodies when tested by FDA accepted methods. No known test method can assure that a product derived from human source materials does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Storage and Stability

Store the control at 2°C–8°C. When stored unopened at 2°C–8°C, the control is stable until the expiration date stated on the label. The control has an open vial stability of one year when stored at 2°C–8°C and one month when stored at room temperature (18°C–25°C).

Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18°C–25°C), at least 15 minutes, depending on remaining volume. Gently invert the control to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the control when not in use.

Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or instrument being used. Technical updates can be found on our website. Safety Data Sheet (SDS) available for professional users at quantimetrix.com.

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Total 25-Hydroxyvitamin D	Units	Level 1 / Lot 169161		Level 2 / Lot 169162		SI Units ¹	Level 1 / Lot 169161		Level 2 / Lot 169162		
		mean	expected range	mean	expected range		mean	expected range	mean	expected range	
Abbott Architect	ng/mL	15.74	12.59 – 18.89	40.81	32.65 – 48.97	nmol/L	38.70	30.96 – 46.44	100.36	80.29 – 120.43	
Beckman Coulter Access Immunoassay System*	ng/mL	31.92	25.53 – 38.30	79.72	63.78 – 95.67	nmol/L	78.49	62.79 – 94.19	196.05	156.84 – 235.26	
Diasorin Liaison*	ng/mL	36.23	28.99 – 43.48	64.66	51.73 – 77.59	nmol/L	89.10	71.28 – 106.92	159.00	127.20 – 190.80	
Euroimmun Analyzer I* ⁴	ng/mL	21.97	17.57 – 26.36	61.93	49.55 – 74.32	nmol/L	54.02	43.21 – 64.82	152.30	121.84 – 182.76	
LC-MS/MS ²		Data not available at this time.									
LC-MS/MS-ID Isotope Dilution* ³		Data not available at this time.									
ROCHE cobas® e 411*	ng/mL	30.34	24.27 – 36.41	71.87	57.50 – 86.25	nmol/L	74.61	59.69 – 89.53	176.74	141.39 – 212.09	
SIEMENS Dimension® ExL	ng/mL	11.99	9.59 – 14.39	38.68	39.95 – 46.42	nmol/L	29.49	23.59 – 35.39	95.13	76.10 – 114.15	

Footnotes

- SI International System of Units
- LC-MS-MS results may vary based on the method used. Values are based on mean results from participating laboratories.
- Measurement procedures are listed by the Joint Committee for Traceability in Laboratory Medicine (JCTLM Method#C8RMP4 and C8RMP3) and traceable to NIST reference material SRM 972.
- 25-OH Vitamin D ELISA (EQ 6411-9601)

* Mean and range based on limited data. Each lab should establish its own mean and range.

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Quantimetrix Corporation
MDSS GmbH