# Dropper® Spinal Fluid Control / Level 1 & 2

**LOT** 45291 45292 **REF** Level 1: 1451-31 Level 2: 1452-31 **2019-05** **CE** **IVD** **!** **2°C-8°C** **i**

**CE** European Conformity / CE-Konformitätsbescheinigung / Conformité aux normes européennes / Conformità europea / Conformidad europea  
**LOT** Lot Number / Beschriftung / Désignation du lot / Numero di lotto / Denominación de lote  
**IVD** Manufactured by / Herstellert / Fabriqué par / Fabricato da / Fabricado por  
For in vitro diagnostic use / In-vitro Diagnostikum / Pour diagnostic in vitro / Per uso diagnostico in vitro / De uso diagnóstico in vitro  
**Cont.** Biological Risk / Biogefährdung / Risque biologique / Rischio biologico / Peligro biológico  
Contents of lot / Inhalt der Packung / Contenu du coffret / Contenuto della confezione / Contenido del estuche  
**REF** Catalog No. / Bestellnr. / N° de catalogue / Catalogo n. / Nº de catálogo  
**!** Caution, See Product Insert / Achtung, Siehe Packungsbeilage / Attention, voir notice d'utilisation / Attenzione, vedere il foglietto illustrativo del prodotto / Atención, consulte el folleto del producto  
**EC REP** Authorized Representative / Bevollmächtigter / Représentant agréé / Rappresentante autorizzato / Representante autorizado  
**i** Temperature limitation / Temperaturbegrenzung / Limites de température / limite di temperatura / límite de temperatura  
**i** Consult instructions for use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni l'uso / Consulte las instrucciones de uso  
**!** Use by (last day of month) / Verwendbar bis (letzter Tag des Monats) / Utilisable jusqu'à (dernier jour du mois indiqué) / Da utilizzare prima del (ultimo giorno del mese) / Estable hasta (último día del mes)

## English

### Intended Use

The Quantimetrix Dropper Spinal Fluid Control is intended as a means of monitoring spinal fluid analyte methods to validate quantitation of patient samples. Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

### Product Description

The Quantimetrix Dropper Spinal Fluid Control is supplied in two levels. It is a ready-to-use liquid, requiring no reconstitution. The Dropper Spinal Fluid Control is prepared in a human protein matrix fortified to target levels with purified chemicals. Preservatives have been added to inhibit microbial growth.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the serum pool have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

### Storage and Stability

Store the control at 2°C–8°C. When stored at 2°C–8°C, the controls are stable until the expiration date stated on the label. Once opened, the controls will be stable for 3 months when stored tightly capped at 2°C–8°C. Discard the controls if turbid or if there is evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

### Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18°C–25°C), at least 15 minutes, depending on remaining volume. Invert gently to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements for the test method. Immediately recap the control and return it to 2°C–8°C when not in use.

### Expected Range

Expected values for the listed lots of controls have been established from interlaboratory data using instrument manufacturers' reagents. Individual laboratory means should fall within the ranges listed. These values should be used as a guide in evaluating the performance of the test methods. Each laboratory should establish its own precision parameters for the methods used to measure each analyte. Mean values and expected ranges apply to all models of the instrument listed unless otherwise noted. A quality control program, Quantrol® is available from Quantimetrix. Please inquire.

### Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Some tests may not have been validated by the manufacturer for use with CSF samples. For information on the instrument specifications, refer to the instrument manufacturer's documentation. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

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- 3 Siemens, Inc., Dimension® EXL, VISTA 500
Quantimetrix Corporation
MDSS GmbH