
Quantimetrix® Dipper^{POCT}® Urinalysis Dipstick Control

Level

Cont. 31 Pouches** - Level 1
31 Pouches** - Level 2

REF 1400-01 **LOT** 241411 2027-02-28

U.S. Pat. No. 9,835,562 U.S. Pat. No. 10,712,283

Quantimetrix
2005 Manhattan Beach Boulevard
Redondo Beach, CA 90278 USA

EC REP MDSS GmbH, Schiffgraben 41
30175 Hannover, Germany

UK Responsible Person: MDSS-UK RP Ltd.
6 Wilmslow Rd., Rusholme,
Manchester M14 5TP
United Kingdom

English

Intended Use

The Quantimetrix Dipper POCT® Urinalysis Dipstick Control is an assayed quality control material used to monitor the performance of urinalysis test strips containing various combinations of reagent pads for bilirubin, blood, creatinine, glucose, ketones, leukocytes, microalbumin, nitrite, pH, protein, specific gravity, and urobilinogen. Dipper POCT may also serve as a control for confirmatory tests such as **K-CHECK**, and **Ictotest**® reagent tablets as well as for **hCG** lateral flow methods.

Summary and Explanation

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

Product Description

The Dipper POCT is a two level urinalysis dipstick control. Each level is contained in a single-use plastic pouch filled with 1.5mL of liquid control material that is ready-to-use, requiring no reconstitution or dilution.

The controls are formulated from a simulated urine matrix fortified with preservatives and stabilizers to maintain product integrity and inhibit microbial growth. Each control level is formulated to target levels with compounds that produce the desired reaction when tested by the methods indicated in the **Intended Use** section.

Caution

Contains human chorionic gonadotropin (hCG) from pregnancy urine. The human hCG source material used in the manufacture of this product has been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C, and HIV 1&2 antibody when tested with FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, 1,2-Propylene Glycol, Level 1; 2,4-Pentanedione, Level 2.

H317 – May cause an allergic skin reaction

P261 – Avoid breathing vapors, mist, or spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The Dipper POCT Urinalysis Dipstick Control kit should be stored at 2°C–8°C and is stable until the expiration date stated on the label.

Do not freeze. Alternatively, each kit may be stored at room temperature (18°C–25°C) for up to 3 months. Kits and/or loose pouches stored at room temperature should be marked with the 3 month expiration date per each laboratory's standard practices. Once opened, pouches should be used within 30 minutes.

Procedure for Dipstick Testing

1 Remove a Dipper POCT pouch from the kit. If the kits have been stored at 2°C–8°C, allow the pouch to equilibrate to room temperature for at least 10 minutes prior to use.

2) Remove a urine dipstick from the container and set aside in an easily accessed location.

3) Hold the Dipper POCT pouch upright and allow the fluid to drain into the lower portion of the pouch. A gentle tap or shake may be required to move the fluid. Become familiarized with the peelable chevron seal and the STOP-PEEL POINT indicated by the Top of the Ribbed Area. (Figure 1.)

4) Use two hands to separate the clear flap from the opaque flap. (Figure 2.)

5) While opening, take care not to apply pressure to the fluid channel as the contents will squirt out. (Figure 3.)

6) Grab each flap placing your thumb on the inside center of the flaps. Apply gentle and consistent pressure to peel open the pouch until you reach the STOP-PEEL POINT. Once you've reached the STOP-PEEL POINT, you will not be able to peel the package further. A triangular shaped opening along the clear flap should be exposed. (Figure 4.)

7) Insert the dipstick and visually verify that all of the reagent pads have been immersed into the fluid. When withdrawing the dipstick, take care to not drag the reagent pads across the clear face of the pouch. The edge of the strip may be drawn along the corner of the pouch to remove excess fluid. Discard the pouch and follow the manufacturer's instructions for reading the dipstick. (Figure 5.)

Procedure for hCG and Confirmatory Tablet Testing

To use the Dipper POCT control for hCG lateral flow and confirmatory tablet methods, begin by following steps 1-6 above. Once the pouch has been opened, a pipette or similar device may be inserted into the pouch to remove the required volume per the hCG test kit or confirmatory tablet test manufacturer's instructions. Discard each pouch after each usage.

Expected Values

All expected value ranges and assignments are determined using pooled data from the system manufacturer and/or internal data generated as per the manufacturer's instructions.

Limitations

Each pouch of Dipper POCT is single use only and should be discarded after each use. Once the control fluid is removed for hCG or confirmatory tablet testing, that pouch should not be used for any further testing. Similarly, once the pouch is used for a urine test strip, it should not be used for any further testing. A small amount of sediment may be present in the pouches. This does not affect the performance of the product.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Detailed information on the limitations of each test method is included in the limitations section of the manufacturers' package insert. Technical updates can be found on our website. The Quality Control Log can be downloaded from the Quantimetrix website at quantimetrix.com or contact Tech Support at (310) 536-0006, option 3.

Siemens Multistix Users

Colors produced by the **urobilinogen** and/or **bilirubin** reactions on these dipsticks with the urinalysis dipstick control may not be characteristic of those shown on the manufacturer's label when reading the dipstick reactions visually. The urobilinogen reactions are consistent and intensify with the increase in the urobilinogen concentration but may not provide an exact color match to those displayed on the label.

Note: Siemens® CLINITEK 50 and Siemens® STATUS or CLINITEK STATUS PLUS may see an Albumin/Creatinine ratio result of "Abnormal" with the Level 1 control.

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


Analytes	Level 1 - Lot 241411	Level 2 - Lot 241412
ROCHE Chemstrip 101 • Urilux S • ROCHE Urisys 1100 (USA) †		
Specific Gravity	1.020 - 1.030	1.005 - 1.020
pH	5.0 - 6.0	7.0 - 9.0
Leukocytes	Negative	25 - 500 Leu/μL (1+ - 3+)
Nitrite	Negative	Positive
Protein	Negative ⁷	30 - 500 mg/dL (1+ - 3+)
Glucose	Normal	100 - ≥1000 mg/dL (1+ - 3+)
Ketones	Negative	15 - 150 mg/dL (1+ - 3+)
Urobilinogen	Normal	1 - 8 mg/dL ^{8*} (1+ - 3+)
Bilirubin	Negative	1 - 6 mg/dL (1+ - 3+)*
Blood	Negative	Trace - 250 Ery/μL (1+ - 3+)
ROCHE cobas u411 (USA) †		
Blood	Negative	50 - 250 Ery/μL (3+ - 5+)
Bilirubin	Negative	1 - 6 mg/dL (1+ - 3+)*
Urobilinogen	Normal	1 - 12 mg/dL (1+ - 4+) ^{8*}
Ketones	Negative	5 - 150 mg/dL (1+ - 4+)
Glucose	Normal	100 - ≥1000 mg/dL (2+ - 4+)
Protein	Negative ⁷	30 - 500 mg/dL (2+ - 4+)
Nitrite	Negative	Positive
Leukocytes	Negative	25 - 500 Leu/μL (1+ - 3+)
pH	5.0 - 6.0	6.5 - 9.0
Specific Gravity	1.020 - 1.030	1.005 - 1.020
Siemens® VISUAL TESTING (Visual Test Strips Only)		
Glucose	Negative	100 - 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Sm - Lg (1+ - 3+)
Ketones	Negative	5 - 160 mg/dL (Sm - Lg)
Specific Gravity	1.020 - 1.030	1.005 - 1.015
Blood	Negative	Small - Lg (1+ - 3+)
pH	5.0 - 6.5	7.5 - 8.5
Protein	Negative	30 - ≥2000 mg/dL (1+ - 4+)
Urobilinogen	Normal (0.2 E.U./dL) ⁷	1.0 - 8.0 E.U./dL*
Nitrite	Negative	Positive
Leukocytes	Negative	Sm - Lg (1+ - 3+)
Creatinine ³	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK 50		
Glucose	Negative	100 - ≥1000 mg/dL (Trace - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - ≥80 mg/dL (Trace - 3+)
Specific Gravity	1.020 - ≥1.030	≤1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5.0 - 6.0	7.5 - ≥9.0
Protein	Negative	30 - ≥300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - 4.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Trace - 3+)
Microalbumin ²	10 - 30 mg/L	30 - 300 mg/L
Creatinine ⁴	10 - 200 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK 500		
Glucose	Negative	100 - 500 mg/dL (Trace - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - ≥80 mg/dL (Trace - 3+)
Specific Gravity	1.020 - ≥1.030	≤1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5.0 - 6.5	7.5 - ≥9.0
Protein	Negative	30 - ≥300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	2.0 - ≥8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Trace - 3+)
Creatinine ³	10 - 100 mg/dL	100 - 300 mg/dL

Analytes	Level 1 - Lot 241411	Level 2 - Lot 241412
Siemens® CLINITEK ADVANTUS		
Glucose	Negative	100 - 500 mg/dL (Trace - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - ≥80 mg/dL (Trace - 3+)
Specific Gravity	1.020 - ≥1.030	≤1.005 - 1.025
Blood	Negative	Small - Large (1+ - 3+)
pH	5.0 - 6.5	7.5 - ≥9.0
Protein	Negative	30 - ≥300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - ≥8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Trace - 3+)
Creatinine ³	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK STATUS or Siemens® CLINITEK STATUS PLUS		
Glucose	Negative	100 - 500 mg/dL (Trace - 2+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	Trace - ≥160 mg/dL (Trace - 4+)
Specific Gravity	1.020 - ≥1.030	1.010 - 1.025
Blood	Negative	Small - Large (1+ - 3+)
pH	5.0 - 6.0	7.5 - ≥9.0
Protein	Negative	30 - ≥300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	2.0 - ≥8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Trace - 3+)
Microalbumin ²	10 - 30 mg/L	30 - 300 mg/L
Creatinine ⁴	10 - 100 mg/dL	100 - 300 mg/dL
hCG	Negative	Positive
Uriscan™ 10 SGL Strips • Visual		
Blood	Negative	10 - 250 RBC/μL (1+ - 3+)
Bilirubin	Negative	0.5 - 3.0 mg/dL (1+ - 3+)
Urobilinogen	Negative (0.1) - Normal (<1)	1 - 12 mg/dL (1+ - 4+)
Ketones	Negative	5 - 100 mg/dL (± - 3+)
Protein	Negative	30 - 1000 mg/dL (1+ - 4+) ⁸
Nitrite	Negative	Positive
Glucose	Negative	200 - 2000 mg/dL (1+ - 4+)
pH	5.0 - 6.0	7.5 - 8.5
Specific Gravity	1.020 - 1.030	1.005 - 1.015
Leukocytes	Negative	25 - 500 WBC/μL (1+ - 3+)
Uriscan™ Optima Urine Analyzers 10 SGL Strips		
Blood	Negative	10 - 250 RBC/μL (1+ - 3+)
Bilirubin	Negative	0.5 - 3.0 mg/dL (1+ - 3+)
Urobilinogen	Normal	1 - 8 mg/dL (1+ - 3+)
Ketones	Negative	5 - 100 mg/dL (± - 3+)
Protein	Negative	30 - 300 mg/dL (1+ - 3+) ⁸
Nitrite	Negative	Positive
Glucose	Negative	100 - 2000 mg/dL (± - 4+)
pH	5.0 - 6.0	7.5 - 8.5
Specific Gravity	1.015 - 1.030	1.000 - 1.015
Leukocytes	Negative	25 - 500 WBC/μL (1+ - 3+)

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Dipper and Dipper POCT are registered trademarks of Quantimetrix.

Quantimetrix Corporation
2005 Manhattan Beach Blvd.
Redondo Beach, CA 90278
+1.310.536.0006

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