Lipasure® Serum Lipoprotein Control / Level 1

LOT 24223
REF 48-7060
2027-01-31
CE
IVD
2°C-8°C
i
Rx only

CE
LOT
Manufactured by Hengstli vers
For in vitro diagnostic use
Biological Risk
Contents of kit
Catalog No.
Caution, See Product Insert
Authorized Representative
Temperature limitation
Consult instructions for use
Use by (last day of month)
U.S. Federal law restricts this device to sale by or only on the order of a licensed healthcare practitioner.

English

Intended Use and Indications for Use

The Quantimetrix Lipasure Lipoprotein Control is intended for monitoring lipoprotein fractions (VLDL to HDL) as well as subfractions (Mid-bands, LDL subfractions 1, 2 and 3) using the Lipoprint LDL System to validate the separation and quantitation of lipoproteins in patient serum samples. Control materials with known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes inter- and intralaboratory parameters for accuracy and precision of the test method.

Product Description

Lipasure is prepared from human serum and supplied lyophilized. It is assayed for VLDL, Midbands C, B and A, LDL subfractions, HDL and Total LDL using the Quantimetrix Lipoprint LDL System. The lyophilized material contains stabilizers and preservatives. The control is provided as a set of 4 x 0.5 mL vials (0.5 mL after reconstitution).

FOR IN VITRO DIAGNOSTIC USE ONLY

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. The Lipasure Control is prepared in a human serum matrix and contains other human source materials. All blood donor units comprising the serum pool have been tested and found nonreactive for Hepatitis B surface antigen, Hepatitis C and HIV 1 & 2 antibodies when tested by FDA accepted methods.

No known test method can assure that a product derived from human source materials does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Storage and Stability

- When stored at 2°C–8°C the unopened, lyophilized material is stable until the expiration date stated on the vial label.
- The reconstituted material is stable for 5 days when kept at 2°C–8°C.
- Discard the control if it is turbid or if it exhibits any signs of microbial contamination.

Materials Needed but Not Provided

- Deionized or distilled water
- A 0.5 mL pipet

Procedure

- Reconstitute the control with 0.5 mL of purified water. Gently swirl to dissolve.
- Allow the reconstituted control to stand for at least 10 minutes at ambient temperature before use.
- Invert gently to assure complete homogeneity of the solution. Avoid foaming.
- Handle the control the same way you would a patient sample. Use as indicated in the Lipoprint LDL System package insert.
- After use immediately recap the control and return it to refrigerated storage (2°C–8°C).
- Include at least one control sample with each run of patient samples.

Qualitative Results

The control will generate a normal lipoprotein profile when used with the Lipoprint LDL System. It will typically exhibit bands for VLDL, Midbands-C, B and A (these comprise IDL), LDL-1, 2 and 3 and HDL when using the Lipoprint LDL template. See Figure 1 below.


Figure 1. Example of Lipoprint band pattern for Lipasure Control using the template

Quantitative Results

The gel tubes are scanned and the relative area of each lipoprotein subfraction is determined from the densitometric scan as described in the Lipoprint LDL System package insert. A representative scan is depicted in Figure 2.


Figure 2. Example of a typical control profile

Expected Values

The assigned values have been established from repetitive runs of the control on the Lipoprint LDL System. See Table on reverse side of insert.

Limitations

- Use of this control with other electrophoresis test systems or precipitation methods for the determination of LDL and HDL may give values other than those assigned for the Lipoprint LDL System.
- Lipasure is not intended for use as a control for total cholesterol or triglycerides.
- Lipasure Level 1 is not intended for use as a control for LDL subfractions 5-7.

Assigned values by Lipoprint LDL System

	Units mg/dL ²		SI Units ¹ mmol/L ³	
	Mean	Range	Mean	Range
VLDL	15	10 – 20	0.394	0.261 – 0.526
Mid-C	14	11 – 17	0.369	0.291 – 0.448
Mid-B	6	2 – 9	0.145	0.053 – 0.236
Mid-A	7	2 – 12	0.180	0.040 – 0.320
LDL-1	23	18 – 27	0.587	0.467 – 0.706
LDL-2	7	0 – 15	0.184	0.000 – 0.388
LDL-3	0	0 – 2	0.002	0.000 – 0.052
LDL-4	0	0 – 2	0.000	0.000 – 0.052
HDL	28	23 – 35	0.717	0.595 – 0.905
Total LDL ⁴	57	45 – 68	1.465	1.172 – 1.757

1 SI International System of Units
 2 Based on total cholesterol of 100 mg/dL

3 Based on total cholesterol of 2.59 mmol/L
 4 Total LDL equals the sum of all Midbands and LDL subfractions present

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Quantimetrix Corporation
MDSS GmbH