# Bilirubin Control – Pediatric / Level 1 & 2

LOT 264041 264042 REF 1331-31 1332-31 2027-06-30 CE IVD 2°C-8°C R<sub>x</sub>only

European Conformity CE-Konformitätserklärung Conformité aux normes européennes Conformité europea Conformidad europea

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U.S. Federal law restricts this device to sale by or only on the order of a licensed healthcare practitioner. Conformité aux normes européennes Conformité europea Conformidad europea

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## English Intended Use

The Quantimetrix Pediatric Bilirubin Control is intended as a means of monitoring various bilirubin assay methods to validate quantitation of patient samples.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. For professional use only.

## Product Description

The Quantimetrix Pediatric Bilirubin Control is supplied as a liquid, 3 x 3 mL vials per box. It is prepared from purified bilirubin in a human protein base. Stabilizers and preservatives have been added.

## Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the protein pool have been tested and found non-reactive for Hepatitis B surface antigen, and for HCV antibody and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is suggested that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

## Storage and Stability

The control should be stored at 2°C–8°C. When stored unopened at 2°C–8°C, the control is stable until the expiration date stated on the label. If the rubber stopper is removed or air is introduced into the vial, the stability of the control is limited to twenty-one (21) days at 2°C–8°C. **Longer stability may be achieved by removing the sample with a syringe.** Discard the control if it is turbid, if there is any evidence of microbial contamination or if it exhibits any green color. Discard controls in the same manner as other biological specimens, according to local guidelines.

## Procedure

Remove the control from the refrigerator and allow it to come to room temperature (18°C–25°C) for 10 to 15 minutes. Gently invert the control to assure homogeneity of the contents. Avoid foaming the control. For maximum stability, lift up the center tab of the metal seal without removing the seal. Pierce the rubber stopper with a syringe and withdraw the required quantity of control. **Do not introduce air into the vial.** Do not return unused material to the vial. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method.

## Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,® is available from Quantimetrix. Please inquire.

## Limitations

Bilirubin is sensitive to light and oxygen. Exposure to these factors may decrease the values from those stated in the expected ranges.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

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Beckman Coulter, Inc., AU Clinical Chemistry System
Roche Diagnostics Corporation, Roche Cobas® 6000, COBAS® INTEGRA 400 Plus
Quantimetrix Corporation
MDSS GmbH