Microbumin[®] Microalbumin Control / Level 1 & 2

LOT 138851
138852 **REF** 1250-01 **2024-05-31** **CE** **IVD** **2°C-8°C** **i**

CE **LOT** **IVD** **Cont.** **REF** **EC REP** **i**

European Conformity
CE-Konformitätsbescheinigung
Conformità aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabbricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnostico in vitro

Biological Risk
Biosigefährdung
Risque biologique
Rischio biologico
Peligro biológico

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English

Intended Use

Quantimetrix's Microbumin Microalbumin Control is intended as a means of monitoring various microalbumin assay methods to validate quantitation of patient samples. Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

Quantimetrix's Microbumin Microalbumin Control is supplied in two levels, 5 x 7 mL each level per box. They are ready to use, liquid, requiring no reconstitution. They are prepared from human urine, fortified to target levels with human albumin and creatinine. Preservatives have been added to inhibit microbial growth.

Warnings and Precautions

POTENTIAL BIOHAZARDOUS MATERIAL. All blood donor units comprising the source plasma used in the manufacture of the albumin have been tested and found non-reactive for Hepatitis B Surface Antigen and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Center for Disease Control's Bio-Safety Level 2 recommendations. Dispose of carefully. Sodium azide may accumulate in plumbing traps and pose a threat of explosion.

Warning Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone.
H317 - May cause an allergic skin reaction.
P261 - Avoid breathing vapors, mist, or spray.
P272 - Contaminated work clothing should not be allowed out of the workplace.
P280 - Wear protective gloves, protective clothing, and eye protection.
P302+P352 - IF ON SKIN: Wash with plenty of water.
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 - Take off contaminated clothing and wash it before reuse.
P501 - Dispose of contents/container in accordance with local, regional, national, and international regulations.
Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The controls should be stored at 2°C–8°C. When stored at 2°C–8°C, the controls are stable until the expiration date stated on the label. When stored at 2°C–8°C between each use, the controls are stable for six months after opening. When using the control with the Micral Test Strips, the control will remain stable for six months after opening, until after 10 uses, or until the expiration date, whichever occurs first. Discard the controls if turbid or if there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18°C–25°C), about 15 - 30 minutes. Invert gently to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the controls and return to 2°C–8°C when not in use.

Expected Values

Expected values for the listed lots of controls have been established from interlaboratory data using instrument manufacturers' reagents. These values should be used as a guide in evaluating the performance of the test methods. Each laboratory should establish its own precision parameters for the methods used to measure each analyte. A quality control program, Quantrol,[™] is available from Quantimetrix. Please inquire.

Mean values and expected ranges apply to all models of the instrument listed unless otherwise noted.

Limitations

The expected mean and ranges were established using instrument manufacturer's reagents available at the time of assay. Any future changes made by the manufacturer of a test method may give different values from those previously recovered. Use of methods other than the ones used to establish the expected values may give different values from the ones indicated. Limitations of the test method are included in the package insert for the reagent or instrument being used.

Depending on the instrument and the reagents used to measure creatinine, the mean creatinine values listed may decrease up to 10% over the entire shelf life of the control.

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Analytes/ Instruments	Method	Units	Level 1 – 138851		Level 2 – 138852	
			mean	expected range	mean	expected range
Microalbumin						
77 Elektronika (Visual/Analyzers)	Reagent Strip	mg/L	≤10	N/A	N/A	30 – 150
Analyticon (CombiScreen Strips Visual)	Reagent Strip	mg/L	10	N/A	N/A	30 – 150
Analyticon (Urilyzer® 100 Pro and CombiScreen Strips)	Reagent Strip	mg/L	10	N/A	N/A	30 – 150
Abbott ARCHITECT	Turbidimetric Immunoassay	mg/L	9.1	7.1 – 11.2	69.6	55.7 – 83.5
Beckman Coulter AU Instruments	Turbidimetric Immunoassay	mg/L	10.5	7.4 – 13.7	68.8	48.2 – 89.5
Roche Chemstrip Micral™	Dipstick Immunoassay	mg/L	N/A	0.0 – 20.0	N/A	50.0 – 100.0
Roche Hitachi/Modular®	Turbidimetric Immunoassay	mg/L	11.0	8.8 – 13.2	64.0	51.2 – 76.8
Roche Hitachi® (Sekisui Rgt)	Turbidimetric Immunoassay	mg/L	8.6	6.7 – 10.5	58.5	45.6 – 71.4
Siemens CLINITEK	Clinitek Microalbumin Rgt Strips	mg/L	N/A	0.0 – ≤10	N/A	30.0 – 80.0
Siemens Dimension®	Turbidimetric Immunoassay	mg/L	10.1	7.0 – 13.1	70.1	56.1 – 84.1
Siemens DCA 2000/Vantage	Turbidimetric Immunoassay	mg/L	12.4	9.9 – 14.9	76.1	60.9 – 91.3
Creatinine						
77 Elektronika (Visual/Analyzers)	Reagent Strip	mg/dL	N/A	10.0 – 100.0	N/A	100.0 – 300.0
Analyticon (Combi Screen Strips Visual)	Reagent Strip	mg/dL	N/A	10.0 – 50.0	N/A	100.0 – 200.0
Analyticon (Urilyzer® 100 Pro and CombiScreen Strips)	Reagent Strip	mg/dL	N/A	10.0 – 50.0	N/A	50.0 – 200.0
Abbott ARCHITECT	Turbidimetric Immunoassay	mg/dL	31.9	25.6 – 38.3	239.2	191.4 – 287.1
Beckman Coulter AU Instruments	Turbidimetric Immunoassay	mg/dL	32.0	22.4 – 41.6	203.0	142.1 – 263.9
Roche Hitachi/Modular®	Turbidimetric Immunoassay	mg/dL	34.2	27.3 – 41.0	236.5	189.2 – 283.8
Roche Hitachi® (Sekisui Rgt)	Picrate	mg/dL	33.5	26.8 – 40.2	230.9	184.7 – 277.1
Roche Hitachi® (Sekisui Rgt)	Enzymatic	mg/dL	36.7	29.3 – 44.0	253.3	202.6 – 303.9
Siemens CLINITEK	Clinitek Microalbumin Rgt Strips	mg/dL	N/A	10.0 – 50.0	N/A	200.0 – 300.0
Siemens Dimension®	Turbidimetric Immunoassay	mg/dL	30.4	24.3 – 36.5	241.7	193.3 – 290.0
Siemens DCA 2000/Vantage	Benedict/Behre	mg/dL	38.9	31.2 – 46.7	261.9	209.5 – 314.3

Footnotes

1 Dimension RXL, EXL
 * Mean and range based on limited data

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
This page intentionally left blank.

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Quantimetrix Corporation
MDSS GmbH