
QuantScopics® Urinalysis Microscopics Control Level 1 (Low) / Level 2 (High)

Insert Revision 09/2021

English Intended use

The QuantScopics Urinalysis Microscopics Control is intended for use as quality control material for the microscopic evaluation of various urine sediment methods.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Microscopic QC controls should be run each day the test is performed.¹

Summary and Explanation

Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.² Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts and crystals.^{3,4}

Product Description

The QuantScopics Urinalysis Microscopics Control is supplied in two levels, 4 x 120 mL bottles per level per box. It is liquid, ready-to-use, requiring no reconstitution or dilution. The control is made from human urine to which stabilized human red blood cells, human white blood cells, and uric acid, calcium oxalate and calcium phosphate crystals have been added. Preservatives have been added to maintain product integrity.

Caution

Contains human urine and human blood cells. All blood donor units comprising the human cell source material used in the manufacture of this product have been tested and found nonreactive for Hepatitis B Surface Antigen and Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

H317 – May cause an allergic skin reaction.

P261 – Avoid breathing vapors, mist, or spray.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P321 – Specific Treatment (see section 4 on the SDS).

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com.

Storage and Stability

The QuantScopics Urinalysis Microscopics Control should be stored tightly capped at 2°C–8°C when not in use. Do not freeze. When stored at 2°C–8°C the controls are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the control from the refrigerator and replace the cap on the control bottle with the pour seal cap included in the control box. Allow the control to come to room temperature (18°C–25°C) for approximately 15–90 minutes depending on the volume remaining in the bottle. Mix the control thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Treat the control as you would a patient sample in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.¹ Immediately close the pour seal cap and store the control at 2°C–8°C when not in use.

Expected Ranges

The expected ranges for each type of formed element found in this lot of QuantScopics Control is listed in the table on page 4. These ranges were determined by assay of multiple bottles of the indicated lot by the methods listed. A 12 mL sample volume was used for each testing method. The samples were centrifuged at 400 RCF (relative centrifugal force) for 5 minutes. After centrifugation, urine sediment was resuspended in either ~0.5 or ~1.0 mL of remaining supernatant according to the respective cell count method procedure. The ranges listed are based on the range of elements observed in 10 high power fields. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters.

Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our website. The Quality Control Log can be downloaded from the Quantimetrix website at quantimetrix.com or contact Tech Support at (310) 536-0006, option 3.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

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- 1 Routine Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, NCCLS Documentation GP 16-A, Approved Guideline, (1995), NCCLS, Wayne PA, 19087.
- 2 Bradley GM, Benson ES, Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed., Philadelphia, PA: Saunders, 1974.
- 3 Bologna, CV, Understanding Laboratory Medicine, St. Louis, MO: CV Mosby Co., 1971.
- 4 Dudas HC, Lab Med, 12:765, 1981.

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Analytes/Method	Level 1		Level 2		Analytes/Method	Level 1		Level 2	
	141681	Units	141682	Units		141681	Units	141682	Units
Red Blood Cells					White Blood Cells				
Atellica UAS 800	3 - 13	Part/µL	30 - 84	Part/µL	Atellica UAS 800	2 - 13	Part/µL	34 - 74	Part/µL
77 Elekronika UriSed Variants / COBIO Variants	0 - 30	Part/µL	10 - 140	Part/µL	77 Elekronika UriSed Variants / COBIO Variants	0 - 15	Part/µL	10 - 80	Part/µL
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	data not available at time of printing				Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	data not available at time of printing			
YD Diagnostics Uriscan PluScope	0 - 20	cell/µL	0 - 100	cell/µL	YD Diagnostics Uriscan PluScope	0 - 20	cell/µL	0 - 60	cell/µL
SYSMEX UF-1000i	data not available at time of printing				SYSMEX UF-1000i	data not available at time of printing			
Beckman Coulter IRIS® Diagnostics iQ200 Series Analyzers	1 - 10	Part/µL	6 - 39	Part/µL	Beckman Coulter IRIS® Diagnostics iQ200 Series Analyzers	2 - 12	Part/µL	18 - 101	Part/µL
Kova® Glassic™ Slide 10 with grids	2 - 13	Part/µL	28 - 106	Part/µL	Kova® Glassic™ Slide 10 with grids	2 - 11	Part/µL	22 - 65	Part/µL
Non-grid slides (~0.5 mL resuspended sediment)	0 - 7	Part/µL	3 - 44	Part/µL	Non-grid slides (~0.5 mL resuspended sediment)	0 - 7	Part/µL	2 - 35	Part/µL
Non-grid slides (~1.0 mL resuspended sediment)	0 - 5	Part/µL	2 - 41	Part/µL	Non-grid slides (~1.0 mL resuspended sediment)	0 - 5	Part/µL	1 - 20	Part/µL
Cen-Slide® System	0 - 7	Part/µL	3 - 42	Part/µL	Cen-Slide® System	0 - 7	Part/µL	1 - 37	Part/µL
Slide & Coverslip (~0.5 mL resuspended sediment)	0 - 6	Part/µL	2 - 28	Part/µL	Slide & Coverslip (~0.5 mL resuspended sediment)	0 - 5	Part/µL	2 - 25	Part/µL
Slide & Coverslip (~1.0 mL resuspended sediment)	0 - 5	Part/µL	2 - 26	Part/µL	Slide & Coverslip (~1.0 mL resuspended sediment)	0 - 4	Part/µL	1 - 17	Part/µL
FisherBrand UriSystem DeciSlide	0 - 6	Part/µL	3 - 56	Part/µL	FisherBrand UriSystem DeciSlide	0 - 6	Part/µL	2 - 54	Part/µL
Roche Cobas u701	0 - 25	Part/µL	23 - 72	Part/µL	Roche Cobas u701	0 - 25	Part/µL	33 - 100	Part/µL

QuantScopics® Urinalysis Microscopics Control

Level 1 (Low) / Level 2 (High)

Casts	Level 1 141681	Level 2 141682
Atellica UAS 800	none	none
77 Elektronika UriSed Variants / COBIO Variants	none	none
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	data not available at time of printing	
YD Diagnostics Uriscan PluScope	none	none
SYSMEX UF-1000i	data not available at time of printing	
Beckman Coulter IRIS® Diagnostics iQ200 Series Analyzers	none	none
Kova® Glassic™ Slide 10 with grids	none	none
Non-grid slides (~0.5 mL resuspended sediment)	none	none
Non-grid slides (~1.0 mL resuspended sediment)	none	none
Cen-Slide® System	none	none
Slide & Coverslip (~0.5 mL resuspended sediment)	none	none
Slide & Coverslip (~1.0 mL resuspended sediment)	none	none
FisherBrand UriSystem DeciSlide	none	none
Roche Cobas u701	NA	NA

Crystals	Level 1 141681	Level 2 141682
Atellica UAS 800	may be present	may be present
77 Elektronika UriSed Variants / COBIO Variants	may be present	may be present
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	data not available at time of printing	
YD Diagnostics Uriscan PluScope	may be present	may be present
SYSMEX UF-1000i	data not available at time of printing	
Beckman Coulter IRIS® Diagnostics iQ200 Series Analyzers	may be present	may be present
Kova® Glassic™ Slide 10 with grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Cen-Slide® System	may be present	may be present
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present
FisherBrand UriSystem DeciSlide	may be present	may be present
Roche Cobas u701	NA	NA

Bacteria	Level 1 141681	Level 2 141682
Atellica UAS 800	may be present	may be present
77 Elektronika UriSed Variants / COBIO Variants	occasionally present	occasionally present
Mindray EH-2030, EH-2050A Plus, EH-2050B Plus, EH-2080B, EH-2080C	data not available at time of printing	
YD Diagnostics Uriscan PluScope	NA	NA
SYSMEX UF-1000i	data not available at time of printing	
Beckman Coulter IRIS® Diagnostics iQ200 Series Analyzers	may be present	may be present
Kova® Glassic™ Slide 10 with grids	may be present	may be present
Non-grid slides (~0.5 mL resuspended sediment)	may be present	may be present
Non-grid slides (~1.0 mL resuspended sediment)	may be present	may be present
Cen-Slide® System	may be present	may be present
Slide & Coverslip (~0.5 mL resuspended sediment)	may be present	may be present
Slide & Coverslip (~1.0 mL resuspended sediment)	may be present	may be present
FisherBrand UriSystem DeciSlide	may be present	may be present
Roche Cobas u701	negative	negative

Quantimetrix Corporation
2005 Manhattan Beach Blvd.
Redondo Beach, CA 90278
+1.310.536.0006

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MDSS GmbH