QuantScopics® Urinalysis Microscopics Control Level 1 (Low) / Level 2 (High)

LOT 232481
232482 **REF** 1481-41
1482-41 **2025-12-31** **CE** **IVD** **2**°C **8**°C **Bio** **i** **Rx only**

CE **LOT** **Manufactured by** **For in vitro diagnostic use** **Biological Risk** **Contents of kit** **Catalog No.** **Caution, See Product Insert** **Autoclave/Reprocessing** **Temperature limitation** **Consult instructions for use** **Use by (best day of month)** **U.S. Federal law restricts this device to sale by or only on the order**

English Intended use

The QuantScopics Urinalysis Microscopics Control is intended for use as quality control material for the microscopic evaluation of various urine sediment methods.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Microscopic QC controls should be run each day the test is performed.¹

Summary and Explanation

Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.² Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts and crystals.^{3,4}

Product Description

The QuantScopics Urinalysis Microscopics Control is supplied in two levels, 4 x 120 mL bottles per level per box. It is liquid, ready-to-use, requiring no reconstitution or dilution. The control is made from human urine to which stabilized human red blood cells, human white blood cells, and uric acid, calcium oxalate and calcium phosphate crystals have been added. Preservatives have been added to maintain product integrity.

Caution

Contains human urine and human blood cells. All blood donor units comprising the human cell source material used in the manufacture of this product have been tested and found nonreactive for Hepatitis B Surface Antigen and Hepatitis C and HIV 1 & 2 antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

H317 – May cause an allergic skin reaction.

P261 – Avoid breathing vapors, mist, or spray.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P321 – Specific Treatment (see section 4 on the SDS).

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com.

Storage and Stability

The QuantScopics Urinalysis Microscopics Control should be stored tightly capped at 2°C–8°C when not in use. Do not freeze. When stored at 2°C–8°C the controls are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the control from the refrigerator and replace the cap on the control bottle with the pour seal cap included in the control box. Allow the control to come to room temperature (18°C–25°C) for approximately 15–90 minutes depending on the volume remaining in the bottle. Mix the control thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Treat the control as you would a patient sample in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.¹ Immediately close the pour seal cap and store the control at 2°C–8°C when not in use.

Expected Ranges

The expected ranges for each type of formed element found in this lot of QuantScopics Control is listed in the table on page 4. These ranges were determined by assay of multiple bottles of the indicated lot by the methods listed. A 12 mL sample volume was used for each testing method. The samples were centrifuged at 400 RCF (relative centrifugal force) for 5 minutes. After centrifugation, urine sediment was resuspended in either ~0.5 or ~1.0 mL of remaining supernatant according to the respective cell count method procedure. The ranges listed are based on the range of elements observed in 10 high power fields. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters.

Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our website. The Quality Control Log can be downloaded from the Quantimetrix website at quantimetrix.com or contact Tech Support at (310) 536-0006, option 3.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

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- 1 Routine Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, NCCLS Documentation GP 16-A, Approved Guideline, (1995), NCCLS, Wayne PA, 19087.
- 2 Bradley GM, Benson ES, Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed., Philadelphia, PA: Saunders, 1974.
- 3 Bologna, CV, Understanding Laboratory Medicine, St. Louis, MO: CV Mosby Co., 1971.
- 4 Dudas HC, Lab Med, 12:765, 1981.

Microscopic Evaluation of Urine Sediment • Mikroskopische Beurteilung von Urinsediment • Evaluation microscopique du sédiment urinaire Valutazione microscopica del sedimento urinario • Evaluación microscópica de la sedimentación presente en la orina

Analytes/Method	Level 1 232481		Level 2 232482		Analytes/Method	Level 1 232481		Level 2 232482	
	Units	Units	Units	Units		Units	Units		
Red Blood Cells					White Blood Cells				
Atellica UAS 800/1500	3 - 12	p/µL	23 - 83	p/µL	Atellica UAS 800	3 - 12	p/µL	28 - 74	p/µL
77 Elektronika UriSed Variants / Analyticon Urilyzer® Cell	0 - 30	p/µL	5 - 135	p/µL	77 Elektronika UriSed Variants / Analyticon Urilyzer® Cell	0 - 15	p/µL	10 - 55	p/µL
Beckman Coulter IRIS® Diagnostics iQ200/DxU Series Analyzers	Not Available	p/µL	Not Available	p/µL	Beckman Coulter IRIS® Diagnostics iQ200/DxU Series Analyzers	Not Available	p/µL	Not Available	p/µL
YD Diagnostics Uriscan Plus	0 - 20	p/µL	0 - 100	p/µL	YD Diagnostics Uriscan Plus	0 - 20	p/µL	0 - 60	p/µL
Roche Cobas u701	0 - 25	p/µL	23 - 69	p/µL	Roche Cobas u701	0 - 25	p/µL	28 - 85	p/µL
Teco UriScope 50 / COBIO Variants	0 - 30	p/µL	48 - 109	p/µL	Teco UriScope 50 / COBIO Variants	0 - 15	p/µL	13 - 68	p/µL
Kova Glassic Slide 10 with grids	2 - 10	p/µL	34 - 107	p/µL	Kova Glassic Slide 10 with grids	2 - 9	p/µL	19 - 69	p/µL
Non-grid slides (~0.5 mL)	0 - 7	p/hpf	4 - 37	p/hpf	Non-grid slides (~0.5 mL)	0 - 7	p/hpf	2 - 22	p/hpf
Non-grid slides (~1.0 mL)	0 - 4	p/hpf	3 - 19	p/hpf	Non-grid slides (~1.0 mL)	0 - 5	p/hpf	1 - 12	p/hpf
Slide & Coverslip (~0.5 mL)	0 - 4	p/hpf	2 - 23	p/hpf	Slide & Coverslip (~0.5 mL)	0 - 4	p/hpf	1 - 13	p/hpf
Slide & Coverslip (~1.0 mL)	0 - 4	p/hpf	0 - 14	p/hpf	Slide & Coverslip (~1.0 mL)	0 - 4	p/hpf	0 - 9	p/hpf
Fisher Desideslide	0 - 4	p/hpf	1 - 18	p/hpf	Fisher Desideslide	0 - 5	p/hpf	1 - 10	p/hpf

QuantScopics® Urinalysis Microscopics Control

Level 1 (Low) / Level 2 (High)

Analytes/Method	Level 1 232481	Level 2 232482
Casts		
Atellica UAS 800	None	None
77 Elektronika UriSed Variants / Analyticon Urilyzer® Cell	May be present	May be present
Beckman Coulter IRIS® Diagnostics iQ200/DxU Series Analyzers	Not Available	Not Available
YD Diagnostics Uriscan Plus	None	None
Roche Cobas u701	None	None
Teco UriScope 50 / COBIO Variants	None	None
Kova Classic Slide 10 with grids	None	None
Non-grid slides (~0.5 mL)	None	None
Non-grid slides (~1.0 mL)	None	None
Slide & Coverslip (~0.5 mL)	None	None
Slide & Coverslip (~1.0 mL)	None	None
Fisher Desideslide	None	None
Crystals		
Atellica UAS 800	None	None
77 Elektronika UriSed Variants / Analyticon Urilyzer® Cell	May be present	May be present
Beckman Coulter IRIS® Diagnostics iQ200/DxU Series Analyzers	Not Available	Not Available
YD Diagnostics Uriscan Plus	None	May be present
Roche Cobas u701	None	None
Teco UriScope 50 / COBIO Variants	None	None
Kova Classic Slide 10 with grids	May be present	May be present
Non-grid slides (~0.5 mL)	May be present	May be present
Non-grid slides (~1.0 mL)	May be present	May be present
Slide & Coverslip (~0.5 mL)	May be present	May be present
Slide & Coverslip (~1.0 mL)	May be present	May be present

Analytes/Method	Level 1 232481	Level 2 232482
Bacteria		
Atellica UAS 800	May be present	May be present
77 Elektronika UriSed Variants / Analyticon Urilyzer® Cell	May be present	May be present
Beckman Coulter IRIS® Diagnostics iQ200/DxU Series Analyzers	Not Available	Not Available
YD Diagnostics Uriscan Plus	None	None
Roche Cobas u701	Negative	Negative
Teco UriScope 50 / COBIO Variants	May be present	May be present
Kova Classic Slide 10 with grids	May be present	May be present
Non-grid slides (~0.5 mL)	May be present	May be present
Non-grid slides (~1.0 mL)	May be present	May be present
Slide & Coverslip (~0.5 mL)	May be present	May be present
Slide & Coverslip (~1.0 mL)	May be present	May be present
Fisher Desideslide	May be present	May be present

Quantimetrix Corporation
2005 Manhattan Beach Blvd.
Redondo Beach, CA 90278
+1.310.536.0006

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MDSS GmbH