# QuantTtest<sup>®</sup> Human Protein Standards / Level 2, 3, 4, & 5

**LOT** 242892 242894 242893 242895 **REF** 3410-02 **2027-05-31** **CE** **IVD** 2°C/8°C **i**

<b>CE</b> European Conformity CE-Konformitätsbescheinigung Conformità aux normes européennes Conformità europea Conformidad europea	<b>LOT</b> Lot Number Designation du lot Numero di lotto Denominación de lote	<b>Manufactured by</b> Hergestellt von Fabbricato per Fabricado por	<b>IVD</b> For in vitro diagnostic use In-vitro Diagnosticum Pour diagnostic in vitro Per uso diagnostico in vitro De uso diagnóstico in vitro	<b>Cont.</b> Biological Risk Biosigefährdung Risque biologique Rischio biologico Peligro biológico	<b>REF</b> Contents of kit Inhalt der Packung Contenu du coffret Contenuto della confezione Contenido del estuche	<b>Caution</b> See Product Insert Achtung, Siehe Packungsbeilage Attention, voir notice d'utilisation Attenzione, vedere il foglietto illustrativo del prodotto Atención, consulte el folleto del producto	<b>EC REP</b> Authorized Representative Bevollmächtigter Représentant agréé Rappresentante autorizzato Representante autorizado	<b>Temperature limitation</b> Temperaturbegrenzungen Limites de température limite di temperatura límite de temperatura	<b>Consult instructions for use</b> Gebrauchsanweisung beachten Consulter les instructions d'utilisation Consultare le istruzioni l'uso Consulte las instrucciones de uso	<b>Use by (last day of month)</b> Verwendbar bis letzter Tag des Monats Utilisable jusqu'à (dernier jour du mois indiqué) Da utilizzare prima del (ultimo giorno del mese) Estable hasta (ultimo día del mes)
--	---	--	---	---	--	---	--	---	---	---

## English

### Intended Use

The QuantTtest Human Protein Standards are intended to be used with the QuantTtest Red (Pyrogallol Red method) Total Protein Assay System. The standards are used to establish standard curves and to calibrate manual and automated methods for cerebrospinal fluid and urine total protein assays.

### Product Description

The QuantTtest Human Protein Standards are available in five protein concentrations: 25, 50, 100 and 200 mg/dL and are packaged in 30 mL dropper bottles. They are liquid, ready-to-use, requiring no reconstitution or dilution and are formulated in an aqueous matrix using purified human albumin and human globulin at a ratio of 70% albumin to 30% globulin. Stabilizers and preservatives, including sodium azide, have been added. The QuantTtest human protein standards are referenced to either SIGMA P5369 or NIST SRM 927c and have been tested by the accepted reference method.<sup>1</sup>

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All source serum units used in the manufacture of purified serum proteins have been tested and found non-reactive for HCV, Hepatitis B surface antigen and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus.

It is recommended that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

DISPOSE OF PROPERLY. Sodium azide may form metal azides in plumbing and pose a threat of explosion.

### Storage and Stability

Store the standards at 2-8 C. When stored at 2-8 C between uses, the standards are stable until the expiration date stated on the label. Discard the standards if turbid or any evidence of microbial contamination is present.

### Procedure

Remove the standards from the refrigerator and allow to come to room temperature (18 – 25 C) for approximately 30 minutes. Gently invert the standards to assure homogeneity of the contents. Avoid foaming. Reagents and other materials for the measurement of total protein are required. Since these vary for each method, refer to the test reagent package insert. For instrument applications, refer to the automated analyzer procedure for CSF and urinary protein to determine which standards are needed for calibration of the instrument. Immediately recap the standards and return to 2 – 8 C when not in use.

### Limitations

Total protein methods show variable absorbance responses to different protein types and concentrations. For further information regarding methodology bias refer to the reagent manufacturer's package insert.<sup>2</sup> Any future changes made by Quantimetrix or instrument manufacturers may give different absorbance values or set points from those previously obtained. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website.

### References

<sup>1</sup>Doumas BT, Bayse DD, Carter RJ et al. A candidate reference method for determination of total protein in serum; I. Development and validation; and II. Test for transferability. Clin Chem 1981;27:1642-1650 and 1651-1654. <sup>2</sup>McElderry LA et al. Six methods for urinary protein compared. Clin Chem 1982;28/2:356-360.

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