QuanTtest® Red Total Protein Assay System

21300
 5210-12
 2020-05

European Conformity
CE-Konformitätskennzeichnung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Bezeichnung
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabrique par
Fabricato da
Fabricado por

For in vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnostico in vitro

Biological Risk
Biogefährdung
Risque biologique
Rischio biologico
Peligro biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Bestellnr.
N° de catalogue
Catalogo n.
N° de catálogo

Caution, See Product Insert
Achtung, Siehe Packungsbeilage
Attention, voir notice d'utilisation
Attenzione, vedere il foglietto illustrativo del prodotto
Atención, consulte el folleto del producto

Authorized Representative
Bevollmächtigter
Représentant agréé
Rappresentante autorizzato
Representante autorizado

Temperature limitation
Temperaturbegrenzung
Limites de température
limite di temperatura
límite de temperatura

Consult instructions for use
Gebrauchsanweisung beachten
Consulter les instructions d'utilisation
Consultare le istruzioni l'uso
Consulte las instrucciones de uso

Use by (last day of month)
Verwendbar bis (letzter Tag des Monats)
Utilisable jusqu'à (dernier jour du mois indiqués)
Da utilizzare prima del (ultimo giorno del mese)
Estable hasta (ultimo día del mes)

English

Intended Use

The QuanTtest Red Total Protein Assay System is intended for the quantitative determination of protein in cerebrospinal fluid and urine for both manual and automated systems.

Summary and Explanation

Determination of total protein in urine and cerebrospinal fluid is valuable in the diagnosis of renal and central nervous system disorders respectively. Urinary protein elevations are commonly seen in the following conditions: strenuous exercise, fever and hypothermia, monoclonal gammopathies, nephritis, nephrosis and diabetic nephropathy, and urinary tract infections. Determination of total protein in cerebrospinal fluid aids in the diagnosis of such conditions as meningitis, encephalitis, poliomyelitis, neurosyphilis, CNS tumors and cerebral hemorrhage.¹

The QuanTtest Red Total Protein Assay System provides a simple method of protein quantitation using microliter samples of urine or cerebrospinal fluid. Total protein is measured by the dye binding method using a complex of Pyrogallol Red and Molybdenum Acid (Mo⁺) described by Fujita et al.² At low pH the dye is red and changes to blue when complexed with protein.

Depending on the equipment and sample size used, a linear response is attained up to 350 mg/dL. The Pyrogallol Red dye binding method can be read at three minutes making the test suitable for automated instrumentation. The sensitivity of the Pyrogallol Red method makes the test appropriate for use with body fluids such as cerebrospinal fluid and urine which have protein concentrations too low for precise measurement by other protein methods.

The biuret method is not sensitive enough for protein determinations in urine. Results with turbidimetric methods may vary depending on the type and concentration of the precipitant, the type of protein, the temperature and the time before reading. Some drugs have been shown to interfere with turbidimetric methods.^{3,4} The QuanTtest Red method has a greater range of linearity than the Coomassie Blue method.

Globulin as well as albumin can be measured with the QuanTtest Red Reagent. Surfactants, copper or iron ions interfere with test results. Samples exhibiting hematuria will give falsely elevated values. This method lessens cuvette staining and can be adapted to automated analyzers.

Principle of the Procedure

The assay is based upon the interaction of protein with the Pyrogallol Red-molybdate complex. Pyrogallol Red is combined with molybdenum acid in the Reagent forming a red complex with a maximum absorbance at 470 nm. Under the conditions of the test, in the presence of protein, a blue-purple color is formed and the maximum absorption of the dye changes from 470 nm to 604 nm. The concentration of the total protein in the sample can be measured by reading the absorbance of the sample at 600 nm. The color response is linear and permits accurate protein determinations.

Reagents

1. QuanTtest Red Reagent is a ready-to-use liquid in 250 mL bottles containing Pyrogallol Red, Sodium Molybdate, Succinic Acid, Sodium Benzoate, Sodium Oxalate, Methanol and Surfactant.
2. Protein Standard Solution Set packaged in 30 mL dropper bottles is supplied in four protein concentrations of 25, 50, 100, and 200 mg/dL as an aqueous solution of 70% human serum albumin and 30% human serum globulin in borate buffered saline. (See accompanying Human Protein Standards package insert.) Sodium azide has been added as a preservative.
3. Liquid Spinal Fluid Controls, Levels 1 and 2, supplied 3 mL per bottle. (See accompanying Spinal Fluid Control package insert). Sodium azide has been added as a preservative.
4. Liquid Human Urine Controls, Levels 1 and 2, supplied 10 mL per bottle. (See accompanying Urine Control package insert). Sodium azide has been added as a preservative.
5. Graph Paper
6. Direction Insert (applicable products listed below)

Product	Catalog No.	Contents
QuanTtest Red Reagent	5210-12	1 x 250 mL
QuanTtest Red Total Protein System	2210-02	200 Test
Reagent		2 x 250 mL
Standards 25, 50, 100, 200 mg/dL		1 x 30 mL each
Dropper® Spinal Fluid Control, Level 1 and 2		1 x 3 mL each
Dropper® Urine Chemistry Control, Level 1 and 2		1 x 10 mL each

Materials Required but Not Provided

1. Pipets capable of delivering 20 µL or 50 µL, and 2.5 mL.
2. Reaction tubes: 13 x 100 mm culture tubes are recommended.
3. Spectrophotometer capable of reading absorbance at 600 nm.
4. Cuvette: glass or disposable plastic.
5. Distilled or deionized water.

Instrumentation

The QuanTtest Red Total Protein Assay System requires the use of a manual or automated spectrophotometer capable of measuring absorbance at 600 nm. Refer to the instrument manufacturer's manual for specific instruction regarding its operation and use. Automated analyzer applications are available upon request from Quantimetrix® Technical Service Department.

For In Vitro Diagnostic Use Only

Warnings and Precautions

The QuanTtest Red Reagent is formulated with acidic methanol. Do not take internally. Do not get in eyes, on skin or on clothing. In case of contact, immediately flush eyes with water and get medical attention. Wash skin with soap and water. Do not pipet by mouth.

The Standards and Controls are formulated from human source material. All blood donor units comprising the serum source material have been tested and found non-reactive for HCV, Hepatitis B surface antigen and HIV antibody when tested by FDA accepted methods. **POTENTIAL BIOHAZARDOUS MATERIAL.** No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

Danger Hazard (H) and Precautionary (P) Statements:

Contains Mixture, Hydrochloric acid, Methyl alcohol, Sodium molybdate dihydrate.

H290 - May be corrosive to metals.

H314 - Causes severe skin burns and eye damage.

P234 - Keep only in original container.

P260 - Do not breathe vapors, mist, or spray.

P280 - Wear protective gloves, protective clothing, and eye protection.

P301+P330+P331 - If swallowed: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 - If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P304+P340 - If inhaled: Remove person to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 - If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 - Immediately call a poison center or doctor.

P501 - Dispose of contents/container in accordance with local, regional, national, territorial, provincial, and international regulations. Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Reconstitution of Reagents

Reagent

The QuanTtest Red Reagent needs no reconstitution or dilution. Remove the Reagent from the refrigerator and aliquot a sufficient quantity of Reagent for the number of tests to be performed. Allow the aliquot to come to room temperature (18°C–25°C) before use.

Standards and Controls

Remove the Standards and Controls from the refrigerator and allow to come to room temperature (18°C–25°C), approximately 30 minutes before use. The Standards and Controls need no reconstitution. Gently invert the Standards and Controls to assure homogeneity of the contents. Avoid foaming.

Storage and Stability

Reagent

1. Store the QuanTtest Red Reagent at 2°C–8°C tightly capped. Protect from light.
2. When stored properly the Reagent is stable until the expiration date stated on the label. Do not use the Reagents past the expiration date stated on each Reagent container label.
3. Discard the Reagent if signs of precipitation or cloudiness develop.

Standards and Controls

1. The Standards and Controls should be stored at 2°C–8°C.
2. When stored at 2°C–8°C between uses, the Standards and Controls are stable until the expiration date stated on the label. Do not use the Standards and Controls past the expiration date stated on each label.
3. Discard the standards or controls if turbid or any evidence of microbial contamination is present.

Specimen Handling

1. Remove the samples to be tested from the refrigerator and allow to come to room temperature (18°C–25°C) before testing.
2. Urine or cerebrospinal fluid samples may be suitable for use without additional sample preparation. No special preservatives are required.
3. Should any blood cells or cellular debris be present in either the urine or spinal fluid specimen, centrifuge the specimen prior to testing.
4. Urinary protein may be determined on a 24-hour collection or on a random specimen.
5. Urine specimens may be stored at 2°C–8°C for up to 24 hours or frozen up to 3 months until assayed. CSF specimens may be stored at 2°C–8°C for several days or frozen for months until assayed.

Procedure

The protein concentration of unknown specimens can be determined by finding the absorbance of the specimen on a standard curve and reading the corresponding protein concentration or can be calculated using a one point calibration (see Calculation of Results below).

QuanTtest® Red Total Protein Assay System

To prepare the standard curve, plot the absorbance obtained for each standard against its respective protein concentration in mg/dL.

If using the calculation method to determine the protein concentration of unknown specimens, use the standard that approximates the concentration of the specimen being tested; e.g., when testing normal CSF specimens measure the 50 mg/dL standard and when testing normal urine specimens measure the 25 mg/dL standard. If the protein concentration of the specimen is higher or lower than these concentrations, use the standard which most closely approximates the protein concentration of the specimen.

Protein can be measured using either a 20 or a 50 µL sample size. The 50 µL sample will give a linear range to 150 mg/dL of protein. To extend the range of linearity of the test method to 350 mg/dL of protein, use a 20 µL sample size (see Limitations section below).

Test Procedure

1. Allow QuanTtest Red Reagent, Standards, Controls, and Specimens to come to room temperature (18°C–25°C).
2. Label test tubes for each sample: Reagent Blank, Standards, Controls, and Patient Specimens to be tested.
3. Pipet 50 µL (20 µL) of each sample into the appropriate tube. Use distilled or deionized water for the Reagent Blank tube.
4. Accurately pipet 2.5 mL of the Reagent to each of the tubes. Mix well by vortexing or agitating the tubes.
5. Allow reaction mixture to stand at room temperature for at least 15 minutes before reading the A600. The color of the reaction mixture is stable for up to 60 minutes.
6. Transfer the Reagent blank to a clean cuvette and adjust the spectrophotometer to read zero at 600 nm.

Note: If the cuvette becomes stained in step 6, dye has adsorbed onto the glass and must be cleaned or erroneously high results will occur. After cleaning, thoroughly rinse with distilled or deionized water and dry. Repeat Step 6.

7. Read each sample against the Reagent blank and record the absorbance values obtained.
8. If preparing a standard curve, plot the absorbance of each standard against its protein concentration on the graph paper provided. Draw the "best fit" straight line through the points. Determine the protein concentration of unknown patient specimens by finding the absorbances on the standard curve and reading the corresponding protein concentrations. If calculating the results, follow the example in the Calculation of Results section below.

Automated Instrumentation Methods

The QuanTtest Red Total Protein Assay System has been evaluated for efficacy when each reaction is read at precisely timed 3 minute intervals, making the method suitable for automated instrumentation. Refer to the specific instrument applications for the Pyrogallol Red method for direction of use of this product. Contact Quantimetrix' Technical Service Department or your instrument manufacturer for these applications.

Calculation of Results

The concentrations of urine or CSF proteins are expressed in mg/dL. The protein concentration of unknown samples can be determined by finding the absorbance of the sample on the standard curve and reading the corresponding protein concentration OR can be calculated as follows:

Calculation of Total Protein

Use the following formula to calculate protein concentrations of the patient samples.

$$\text{Protein Concentration (mg/dL)} = \frac{\text{absorbance of sample X Concentration of std. (mg/dL)}}{\text{absorbance of standard}}$$

The following are examples for illustrative purposes only.

Cerebrospinal Fluid

$$\begin{aligned} \text{CSF Patient Specimen } A_{600} &= 0.200 \\ 50 \text{ mg/dL Standard } A_{600} &= 0.310 \\ (0.200 / 0.310) \times 50 \text{ mg/dL} &= 32.3 \text{ mg/dL} \end{aligned}$$

Urine

$$\begin{aligned} \text{Urine Patient Specimen } A_{600} &= 0.120 \\ 25 \text{ mg/dL Standard } A_{600} &= 0.165 \\ (0.120 / 0.165) \times 25 \text{ mg/dL} &= 18.2 \text{ mg/dL} \end{aligned}$$

Calculation of 24 Hour Urinary Protein

1. Measure the 24 hour urine volume (mL).
2. Divide this volume by 100 to determine the volume expressed in deciliters (dL).
3. Multiply the 24 hour urine volume in dL by the protein concentration in mg/dL as determined by QuanTtest Red procedure.
4. The resulting number is the amount of protein expressed as mg per 24 hours.

The following is an example for illustrative purposes only.

$$\begin{aligned} \text{Patient Urine Volume} &= 1,550 \text{ mL} \\ \text{Patient Urinary Protein} &= 18.2 \text{ mg/dL} \\ (1550 \text{ mL} / 100 \text{ mL/dL}) \times 18.2 \text{ mg/dL} &= 282.1 \text{ mg/24 hours} \end{aligned}$$

Quality Control

A quality control program is recommended for all clinical laboratories. The analysis of control material with assayed values in both the normal and abnormal ranges is recommended for monitoring the performance of the procedure. The values recovered for the controls should fall within the laboratory's acceptable range. A Quality Control Program, Quantrol,® is available from Quantimetrix. Please inquire.

Limitations

Linearity of the test method is 150 mg/dL of total protein when using a 50 µL sample size. Linearity may be extended by decreasing the specimen volume to 20 µL. The total protein concentration then becomes linear to 350 mg/dL. Greater accuracy may be obtained by diluting specimens of higher

protein levels to approximately 100 mg/dL with distilled or deionized water, reassaying and correcting for the dilution factor when calculating the total protein concentration.

Some types of surface active agents may affect the color. Cationic surfactants in general, give the same color as proteins. Because anions inhibit the color reaction, wash the equipment thoroughly, using distilled water, until no surface active agent remains. Dry equipment thoroughly before using.

Hemoglobin gives about one-half the color reaction of albumin. Samples exhibiting hematuria will give a falsely elevated value.

Copper and iron ions can cause errors in measurement. Special care should be taken to avoid contamination of the apparatus.

The Reagent contains light sensitive components. Protect from light.

The QuanTtest Red system may underestimate low molecular weight proteins and immunoglobulin lightchains (Bence Jones Proteins).

In routine use, cuvettes, culture tubes and equipment are rarely stained. However, when many specimens are run at the same time or if a specimen containing a high concentration of protein is tested, some of the dye may remain. If this occurs, clean glassware and equipment thoroughly.

Technical updates can be found on our website.

Expected Values

Normal ranges may vary according to patient population and laboratory conditions. Each laboratory should establish its own normal range. Using the Pyrogallol Red method on specimens from clinically healthy subjects, protein levels in cerebrospinal fluid can be expected to range from 15 to 45 mg/dL and up to 60 mg/dL in healthy elderly subjects. Urinary excretion of protein is normally less than 200 mg/24 hours in volumes of 1,000 to 1,500 mL/24 hours. Random normal urine specimens may range from 0 to 20 mg/dL of protein.

Specific Performance Characteristics Specificity

The specificity of the QuanTtest Red test method was determined by analysis of urine samples following the addition of potential interfering substances. Hemoglobin, iron and copper ions were found to interfere. See the table of interfering substances following.

Table of Interfering Substances

Additive Name	Added Conc. (mg/dL)	Measured Conc. of Total Protein (mg/dL)	Amount Normally Present in Urine
None	—	29	—
Urea Nitrogen	2,300	29	10.5 g/day
Uric Acid	300	30.5	250-750 mg/day
Creatinine	300	28.5	male 1000-1700 mg/day female 1175-1468 mg/day
Ammonia	250	30	0.41 - 0.82 g/day
Glucose	5,000	28	30 - 130 mg/day
Ascorbic Acid	200	28.5	0.1 - 0.4 mg/kg/day
Glutathione	50	29	—
Citric Acid	200	29.5	0.3 - 0.9 mg/day
Oxalic Acid	100	26	15 - 30 mg/day
Tartaric Acid	200	28	—
Na	600	30	6 - 8.4 g/day
K	250	28	1.8 - 2 g/day
Cl	1,000	30	11.1 - 18.2 g/day
Ca	100	30	240 - 320 mg/day
Mg	30	29	70 - 155 mg/day
P	200	33	—
Cu	10	39	250 µg/day
Fe ²⁺	10	36.5	60 - 100 µg/day
Bilirubin	20	30.5	3.08 ± 0.28 mg/day

Precision

The precision of the QuanTtest Red System was established by analyzing a series of Quantimetrix Urine and Spinal Fluid Controls. The assay was calibrated with Quantimetrix standards as included in the kit. The reaction was read at 3 minutes and after 15 minutes on a Gilford Star II spectrophotometer.

When read at 3 minutes:


	Within Run (n=10)			Between Run (n=18)		
	Mean (mg/dL)	SD	%CV	Mean (mg/dL)	SD	%CV
Urine Level 1	15.3	2.3	14.9	14.9	1.3	8.7
Urine Level 2	62.5	2.5	4.0	61.8	2.5	4.1
CSF Level 1	36.3	2.9	7.9	35.3	3.3	9.4
CSF Level 2	66.7	2.1	3.2	65.2	3.3	5.1

When read after 15 minutes:

	Within Run (n=10)			Between Run (n=18)		
	Mean (mg/dL)	SD	%CV	Mean (mg/dL)	SD	%CV
Urine Level 1	15.2	0.4	2.6	15.5	1.3	7.5
Urine Level 2	63.3	3.0	4.8	62.8	1.9	3.1
CSF Level 1	36.7	2.6	7.2	35.9	3.4	9.5
CSF Level 2	69.6	2.6	3.7	68.0	2.6	3.9

QuanTtest® Red Total Protein Assay System

Example Graph

Correlation

Comparison of the QuanTtest Red method with a Coomassie Blue method (Quantimetrix) for protein determination gives the following correlation when tested using 29 patient urine samples.

n=29	Mean	Correlation Coefficient (r)	Slope	Intercept
Coomassie Blue	43.04	0.9923	1.05	-0.18
QuanTtest Red	45.20			

References

1. Pesce AJ, Kaplan LA. Methods in Clinical Chemistry. CV Mosby Co, 1987.
2. Fujita Y, Mori I, Kitano S. Bunseki Kagaku 1983;32:379.
3. Lievens MM, Celis PJ. Drug interference in turbidimetry and colorimetry of proteins in urine. Clin Chem 1982;28:2328.
4. Muir A, Hensley WJ. Pseudoproteinuria due to penicillins in the turbidimetric measurement of proteins with trichloroacetic acid. Clin Chem 1979;25:1662-1663.

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Quantimetrix Corporation
MDSS GmbH