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English Intended Use

The quanTtest red Total Protein Assay System is intended for the quantitative determination of protein in cerebrospinal fluid and urine for both manual and automated systems.

Summary and Explanation

Determination of total protein in urine and cerebrospinal fluid is valuable in the diagnosis of renal and central nervous system disorders respectively. Urinary protein elevations are commonly seen in the following conditions: strenuous exercise, fever and hypothermia, monoclonal gammopathies, nephritis, nephrosis and diabetic nephropathy, and urinary tract infections. Determination of total protein in cerebrospinal fluid aids in the diagnosis of such conditions as meningitis, encephalitis, poliomyelitis, neurosyphilis, CNS tumors and cerebral hemorrhage.¹

The quanTtest red Total Protein Assay System provides a simple method of protein quantitation using microliter samples of urine or cerebrospinal fluid. Total Protein is measured by the dye binding method using a complex of pyrogallol red and molybdenum acid (Mo+⁶) described by Fujita et al². At low pH the dye is red and changes to blue when complexed with protein.

Depending on the equipment and sample size used, a linear response is attained up to 350 mg/dL. The Pyrogallol Red dye binding method can be read at three minutes making the test suitable for automated instrumentation. The sensitivity of the Pyrogallol Red method makes the test appropriate for use with body fluids such as cerebrospinal fluid and urine which have protein concentrations too low for precise measurement by other protein methods.

The biuret method is not sensitive enough for protein determinations in urine. Results with turbidimetric methods may vary depending on the type and concentration of the precipitant, the type of protein, the temperature and the time before reading. Some drugs have been shown to interfere with turbidimetric methods.^{3,4} The quanTtest red method has a greater range of linearity than the Coomassie Blue method.

Globulin as well as albumin can be measured with the quanTtest red Reagent. Surfactants, copper or iron ions interfere with test results. Samples exhibiting hematuria will give falsely elevated values. This method lessens cuvette staining and can be adapted to automated analyzers.

Principle of the Procedure

The assay is based upon the interaction of protein with the Pyrogallol Red-molybdate complex. Pyrogallol Red is combined with molybdenum acid in the Reagent forming a red complex with a maximum absorbance at 470 nm. Under the conditions of the test, in the presence of protein, a blue-purple color is formed and the maximum absorption of the dye changes from 470 nm to 604 nm. The concentration of the total protein in the sample can be measured by reading the absorbance of the sample at 600 nm. The color response is linear and permits accurate protein determinations.

Reagents

1. QuanTtest red Reagent is a ready-to-use liquid in 250 mL bottles containing Pyrogallol Red, Sodium Molybdate, Succinic Acid, Sodium Benzoate, Sodium Oxalate, Methanol and Surfactant.
2. Protein Standard Solution Set packaged in 30 mL dropper bottles is supplied in four protein concentrations of 25, 50, 100, and 200 mg/dL as an aqueous solution of 70% human serum albumin and 30% human serum globulin in borate buffered saline. (See accompanying Human Protein Standards package insert). Sodium azide has been added as a preservative.
3. Liquid Spinal Fluid Controls, Levels 1 and 2, supplied 3 mL per bottle. (See accompanying Spinal Fluid Control package insert). Sodium azide has been added as a preservative.
4. Liquid Human Urine Controls, Levels 1 and 2, supplied 10 mL per bottle. (See accompanying Urine Control package insert). Sodium azide has been added as a preservative.
5. Graph Paper
6. Direction Insert (applicable products listed below)

Product	Catalog No.	Contents
QuanTtest Red Reagent	5210-12	1 x 250 mL
QuanTtest Red Total Protein System Reagent	2210-02	200 Test 2 x 250 mL
Standards 25, 50, 100, 200 mg/dL		1 x 30 mL each
Dropper® Spinal Fluid Control, Level 1 and 2		1 x 3 mL each
Dropper® Urine Chemistry Control, Level 1 and 2		1 x 10 mL each

Materials Required but Not Provided

1. Pipets capable of delivering 20 µL or 50 µL, and 2.5 mL.
2. Reaction tubes: 13 x 100 mm culture tubes are recommended.
3. Spectrophotometer capable of reading absorbance at 600 nm.
4. Cuvette: glass or disposable plastic.
5. Distilled or deionized water.

Instrumentation

The quanTtest red Total Protein Assay System requires the use of a manual or automated spectrophotometer capable of measuring absorbance at 600 nm. Refer to the instrument manufacturer's manual for specific instruction regarding its operation and use. Automated analyzer applications are available upon request from Quantimetrix' Technical Service Department.

For In Vitro Diagnostic Use Only

Warnings and Precautions

The quanTtest red Reagent is formulated with acidic methanol. Do not take internally. Do not get in eyes, on skin or on clothing. In case of contact, immediately flush eyes with water and get medical attention. Wash skin with soap and water. Do not pipet by mouth.

The Standards and Controls are formulated from human source material. All blood donor units comprising the serum source material have been tested and found non-reactive for HCV, Hepatitis B surface antigen and HIV antibody when tested by FDA accepted methods. **POTENTIAL BIOHAZARDOUS MATERIAL.** No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

Reconstitution of Reagents

Reagent

The quanTtest red Reagent needs no reconstitution or dilution. Remove the Reagent from the refrigerator and aliquot a sufficient quantity of Reagent for the number of tests to be performed. Allow the aliquot to come to room temperature (18-25° C) before use.

Standards and Controls

Remove the Standards and Controls from the refrigerator and allow to come to room temperature (18-25° C), approximately 30 minutes before use. The Standards and Controls need no reconstitution. Gently invert the Standards and Controls to assure homogeneity of the contents. Avoid foaming.

Storage and Stability

Reagent

1. Store the quanTtest red Reagent at 2-8° C tightly capped. Protect from light.
2. When stored properly the Reagent is stable until the expiration date stated on the label. Do not use the Reagents past the expiration date stated on each Reagent container label.
3. Discard the Reagent if signs of precipitation or cloudiness develop.

Standards and Controls

1. The Standards and Controls should be stored at 2-8° C.
2. When stored at 2-8° C between uses, the Standards and Controls are stable until the expiration date stated on the label. Do not use the Standards and Controls past the expiration date stated on each label.
3. Discard the standards or controls if turbid or any evidence of microbial contamination is present.

Specimen Handling

1. Remove the samples to be tested from the refrigerator and allow to come to room temperature (18-25° C) before testing.
2. Urine or cerebrospinal fluid samples may be suitable for use without additional sample preparation. No special preservatives are required.
3. Should any blood cells or cellular debris be present in either the urine or spinal fluid specimen, centrifuge the specimen prior to testing.
4. Urinary protein may be determined on a 24-hour collection or on a random specimen.
5. Urine specimens may be stored at 2-8° C for up to 24 hours or frozen up to 3 months until assayed. CSF specimens may be stored at 2-8° C for several days or frozen for months until assayed.

Procedure

The protein concentration of unknown specimens can be determined by finding the absorbance of the specimen on a standard curve and reading the corresponding protein concentration or can be calculated using a one point calibration (see Calculation of Results below).

To prepare the standard curve, plot the absorbance obtained for each standard against its respective protein concentration in mg/dL.


If using the calculation method to determine the protein concentration of unknown specimens, use the standard that approximates the concentration of the specimen being tested; e.g., when testing normal CSF specimens measure the 50 mg/dL standard and when testing normal urine specimens measure the 25 mg/dL standard. If the protein concentration of the specimen is higher or lower than these concentrations, use the standard which most closely approximates the protein concentration of the specimen. Protein can be measured using either a 20 or a 50 µL sample size. The 50 µL sample will give a linear range to 150 mg/dL of protein. To extend the range of linearity of the test method to 350 mg/dL of protein, use a 20 µL sample size (see Limitations section below).

Test Procedure

1. Allow quanTtest red Reagent, Standards, Controls, and Specimens to come to room temperature (18-25° C).
2. Label test tubes for each sample: Reagent Blank, Standards, Controls, and Patient Specimens to be tested.

QuantTest Red Total Protein Assay System

Example Graph

Correlation

Comparison of the quanTTest red method with a Coomassie Blue method (Quantimetrix) for protein determination gives the following correlation when tested using 29 patient urine samples.

n=29	Mean	Correlation Coefficient (r)	Slope	Intercept
Coomassie Blue	43.04	0.9923	1.05	-0.18
QuantTest red	45.20			

References

1. Pesce AJ, Kaplan LA. Methods in Clinical Chemistry. CV Mosby Co, 1987.
2. Fujita Y, Mori I, Kitano S. Bunseki Kagaku 1983;32:379.
3. Lievens MM, Celis PJ. Drug interference in turbidimetry and colorimetry of proteins in urine. Clin Chem 1982;28:2328.
4. Muir A, Hensley WJ. Pseudoproteinuria due to penicillins in the turbidimetric measurement of proteins with trichloroacetic acid. Clin Chem 1979;25:1662-1663.
5. Young GS et al. Effect of drugs on clinical laboratory tests. Clin Chem 1975;21:1D.
6. Martin EW. Hazard of Medication (Alexander SF, Farage, DJ, Hassan WE Jr. eds.) Philadelphia PA and Toronto Canada: JB Lippincott Co 1971:169-189.
7. Constantino NV, Kabat HF. Drug-induced modification for laboratory test values, revised 1973. Am J Hosp Pharm 1973;30:24.
8. Watanabe N, Makino K, Kamei S, Okubo A, Yamanaka M, Osawa S. The Japanese Journal of Clinical Pathology 1984;32 suppl:227.
9. Yoshizaki H, Osawa S, Furuya S. The Japanese Journal of Clinical Pathology 1984;32 suppl:227.
10. Bradford MM. A rapid and sensitive method for the quantization of microgram quantities of protein utilizing the principle of protein dye-binding. Anal Biochem 1976;72:248.
11. Kingsbury FB, Clark CP, Williams G, Post AL. J Lab Clin Med 1926;11:981.
12. Meulemans O. Clin Chem, Acta 1960;5:757.
13. Saito M, Kitamura M, Niwa M. Rinsho-kagaku-bunseki II. Tokyo Kagakudojin 1979:21.

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Références

- Pesce AJ, Kaplan LA. Methods in Clinical Chemistry. CV Mosby Co, 1987.
- Fujita Y, Mori I, Kitano S. Bunseki Kagaku 1983;32:379.
- Lievens MM, Celis PJ. Drug interference in turbidimetry and colorimetry of proteins in urine. Clin Chem 1982;28:2328.
- Muir A, Hensley WJ. Pseudoproteinuria due to penicillins in the turbidimetric measurement of proteins with trichloroacetic acid. Clin Chem 1979;25:1662-1663.
- Young GS et al. Effect of drugs on clinical laboratory tests. Clin Chem 1975;21:1D.
- Martin EW. Hazard of Medication (Alexander SF, Farage, DJ, Hassan WE Jr. eds.) Philadelphia PA and Toronto Canada: JB Lippincott Co 1971:169-189.
- Constantino NV, Kabat HF. Drug-induced modification for laboratory test values, revised 1973. Am J Hosp Pharm 1973;30:24.
- Watanabe N, Makino K, Kamei S, Okubo A, Yamanaka M, Osawa S. The Japanese Journal of Clinical Pathology 1984;32 suppl:227.
- Yoshizaki H, Osawa S, Furuya S. The Japanese Journal of Clinical Pathology 1984;32 suppl:227.
- Bradford MM. A rapid and sensitive method for the quantization of microgram quantities of protein utilizing the principle of protein dye-binding. Anal Biochem 1976;72:248.
- Kingsbury FB, Clark CP, Williams G, Post AL. J Lab Clin Med 1926;11:981.
- Meulemans O. Clin Chem, Acta 1960;5:757.
- Saito M, Kitamura M, Niwa M. Rinsho-kagaku-bunseki II. Tokyo Kagakudojin 1979:21.

QuanTtest Red Total Protein Assay System

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