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English

Intended Use

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are assayed controls intended as a means of monitoring the measurement of conductivity and osmolality by cystic fibrosis screening procedures and to validate the analysis of Chloride, Potassium, and Sodium by confirmatory methods used to quantitate patient samples.¹

Product Description

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are liquid, ready to use, requiring no reconstitution or dilution. They are supplied in three levels, nine vials per box, 3 x 3 mL each level. They are prepared in an aqueous simulated human sweat matrix² to which preservatives have been added.

Warnings and Precautions

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing do not contain human source material. No special precautions are required in handling the product other than those routinely used in the laboratory.

Warning Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone. Urea, Sodium Chloride, Potassium Chloride.

H317 – May cause an allergic skin reaction

P261 – Avoid breathing vapors, mist, or spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The Sweat Controls should be stored refrigerated (2°C– 8°C) or at room temperature (18°C– 25°C). Do not freeze. When stored unopened at 2°C– 8°C or at 18°C– 25°C the controls are stable until the expiration date stated on the label. When stored at 2°C– 8°C between each use the controls are stable for 3 months after opening. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the controls from the refrigerator and allow them to come to room temperature (18°C– 25°C), at least 15 minutes, depending on remaining volume. Invert gently to ensure homogeneity of the contents. Immediately recap the control when not in use.

Expected Values

The expected values have been established in the Quantimetrix laboratory and from interlaboratory data using the listed manufacturer's instruments. Use of control materials having known component concentrations are an integral part of diagnostic procedures. Monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Each laboratory should establish its own quality control standards for the test method used. A quality control program, Quantrol, is available from Quantimetrix. Please inquire.

Limitations

Any future changes made by the manufacturer of a test method may give a different value from the one indicated. Limitations of the test method are included in the information provided by the instrument manufacturer. Technical updates can be found on our website.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone. Urea, Sodium Chloride, Potassium Chloride.

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Quantimetrix Corporation
MDSS GmbH