

# Technical Update: Dipper® Urinalysis Dipstick Control

Lot #44551

Tech number: 290-18-44-1

Users: All Dipper Lot #44551 Urinalysis Dipstick Control Users • Date issued: 10/17/18



## English

Please be advised of the following revision to the Dipper® Urinalysis Dipstick Control package inserts: (Table 1)

## German

Bitte nehmen Sie folgende Änderung in der Packungsbeilage für Dipper® Urinalysis Dipstick Control zur Kenntnis: (Tabelle 1)

## French

Veuillez prendre note des correctifs suivants concernant les notices d'utilisation du Dipper® Urinalysis Dipstick Control : (Tableau 1)

## Italian

Si prega di notare la seguente revisione per l'inserito informativo della confezione del Dipper® Urinalysis Dipstick Control: (Tabella 1)

## Spanish

Por favor, tenga en cuenta las siguientes modificaciones en el prospecto del Dipper® Urinalysis Dipstick Control: (Tabla 1)

**Table 1 / Tabelle 1 / Tableau 1 / Tabella 1 / Tabla 1**

	Level 1	Level 2
<b>Confirmatory and Other Tests</b>		
<b>Refractometer (Specific Gravity)</b>	<b>1.019 – 1.025</b>	1.008 - 1.016