

Technical Update: Dropper® Urinalysis Dipstick Control

Lot #44651

Tech number: 043-19-44-1

All Roche Chemstrip UA Analyzers, Roche Criterion, Roche Cobas u411, Roche Urisys 1800 (US) and Roche Midtron M, Roche Midtron Junior and Urisys 1800 (International) of the Dropper Urinalysis Dipstick Control Lot #44651 • Date issued: 2/13/19



English

Please be advised of the following revision to the Dropper® Urinalysis Dipstick Control package inserts: (Table 1)

Deutsch

Bitte nehmen Sie folgende Änderung in der Packungsbeilage für Dropper® Urinalysis Dipstick Control zur Kenntnis: (Tabelle 1)

Français

Veuillez prendre note des correctifs suivants concernant les notices d'utilisation du Dropper® Urinalysis Dipstick Control : (Tableau 1)

Italiano

Si prega di notare la seguente revisione per l'inserito informativo della confezione del Dropper® Urinalysis Dipstick Control: (Tabella 1)

Español

Por favor, tenga en cuenta las siguientes modificaciones en el prospecto del Dropper® Urinalysis Dipstick Control: (Tabla 1)

Table 1 / Tabelle 1 / Tableau 1 / Tabella 1 / Tabla 1

| Analytes | Level 1 - 44651 | Level 2 - 44652 |
|------------------------------------------|-----------------|-----------------|
| ROCHE Chemstrip UA Urine Analyzer | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |
| ROCHE Criterion | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |
| ROCHE cobas u 411 | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |
| ROCHE Urisys 1800 | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |

INTERNATIONAL USE ONLY

| Analytes | Level 1 - 44651 | Level 2 - 44652 |
|-----------------------------------------------|-----------------|-----------------|
| ROCHE Midtron M | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |
| ROCHE Midtron Junior | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |
| ROCHE Urisys 1800 or Roche cobas u 411 | | |
| Specific Gravity | 1.015 - 1.025 | 1.000 - 1.020 |