

# Technical Update: Dropper® Urinalysis Dipstick Control

Lots #44652 / 44722

Tech number: 009-20-44-1

All Siemens® CLINITEK users of Dropper Urinalysis Dipstick Control Lots #44652 / 44722 • Date issued: 01/22/2020



## English

Please be advised of the following revision to the Dropper® Urinalysis Dipstick Control package inserts: (Table 1)

## Deutsch

Bitte nehmen Sie folgende Änderung in der Packungsbeilage für Dropper® Urinalysis Dipstick Control zur Kenntnis: (Tabelle 1)

## Français

Veuillez prendre note des correctifs suivants concernant les notices d'utilisation du Dropper® Urinalysis Dipstick Control : (Tableau 1)

## Italiano

Si prega di notare la seguente revisione per l'inserito informativo della confezione del Dropper® Urinalysis Dipstick Control: (Tabella 1)

## Español

Por favor, tenga en cuenta las siguientes modificaciones en el prospecto del Dropper® Urinalysis Dipstick Control: (Tabla 1)

**Table 1 / Tabelle 1 / Tableau 1 / Tabella 1 / Tabla 1**

Analytes	Level 2 - 44652 / 44722
<b>Siemens CLINITEK 50</b>	
Ketones	Trace - $\geq 80$ mg/dL (Tr - 3+)
<b>Siemens CLINITEK 500</b>	
Ketones	Trace - $\geq 80$ mg/dL (Tr - 3+)
<b>Siemens CLINITEK ADVANTUS</b>	
Ketones	Trace - $\geq 80$ mg/dL (Tr - 3+)
<b>Siemens CLINITEK STATUS</b>	
Ketones	Trace - $\geq 160$ mg/dL (Tr - 4+)