

Technical Update: Microbumin® Microalbumin Control

Lot #25370

Tech number: 032-19-25-1

All Beckman Coulter AU Instrument Users of Microbumin Microalbumin Control Lot #25370 • Date issued: 2/1/19



English

Please be advised of the following revision to the Microbumin® Microalbumin Control package inserts: (Table 1)

Deutsch

Bitte nehmen Sie folgende Änderung in der Packungsbeilage für Microbumin® Microalbumin Control zur Kenntnis: (Tabelle 1)

Français

Veillez prendre note des correctifs suivants concernant les notices d'utilisation du Microbumin® Microalbumin Control : (Tableau 1)

Italiano

Si prega di notare la seguente revisione per l'inserito informativo della confezione del Microbumin® Microalbumin Control: (Tabella 1)

Español

Por favor, tenga en cuenta las siguientes modificaciones en el prospecto del Microbumin® Microalbumin Control: (Tabla 1)

Table 1 / Tabelle 1 / Tableau 1 / Tabella 1 / Tabla 1

Analytes/Instruments	Method	Units	Level 1 – 25371		Level 2 – 25372	
			mean	expected range	mean	expected range
Creatinine						
Beckman Coulter AU Instruments	Turbidimetric Immunoassay	mg/dL	27.4	19.2 – 35.6	194.3	136.0 – 252.6

Footnotes for values	Fußnoten für werte	Apostilles pour des valeurs	Note a piè di pagina per i valori	Notas al pie de la página para los valores
ENGLISH 9 Atypical color	DEUTSCH 9 Atypische Farbe	FRANÇAIS 9 Couleur atypique	ITALIANO 9 Colore atipico	ESPAÑOL 9 Color anormal